

Section Contents

SPE Column Selection	1-002	SPE Accessories	1-051
SPE Procedure	1-005	eVol XR Dispensing System	1-062
SOLA SPE Cartridges and Plates	1-007	QuEChERS Products	1-065
SPE Cartridges and Plates	1-012	Syringe Filters	1-072
Automated SPE	1-048		

Featured Products

SOLA SPE

Maximizing reproducibility and throughput for bio-analysis

➤ **PAGE 1-007**



Syringe Filters

Reliable elimination of particles and micro organisms

➤ **PAGE 1-072**



QuEChERS – Dispersive SPE

Efficient pesticide sample preparation from food matrices

➤ **PAGE 1-065**



SLE

Fast and effective sample clean up of biological samples

➤ **PAGE 1-035**



HyperSep SPE

Reproducible and clean sample extraction

➤ **PAGE 1-012**



Protein Precipitation

Fast removal of proteins from biological samples

➤ **PAGE 1-034**



Thermo Scientific SPE Products

SPE Cartridge Selection

Thermo Scientific™ SPE products have been developed to meet the requirements of today's sample preparation challenges. Thermo Scientific SPE cartridges and 96-well plates are offered in a range of phases and bed weights, ideal for use in application areas such as: Pharmaceutical, Biochemical, Environmental, and Food and Beverage.

The following flow chart is designed to help in selecting the correct phase for the compound of interest.

SPE Phase Selection

Sample Matrix	Solubility of Sample	Polarity of Sample	Separation Mode	Recommended HyperSep Sorbent(s)
Aqueous	Water Soluble	Non-Polar	Reversed Phase	SOLA HRP, C18, C8, Phenyl, Retain PEP
		Moderately Polar	Reversed Phase	SOLA HRP, C18, C8, Retain PEP
		Polar	Reversed Phase	SOLA HRP, Hypercarb Retain PEP
		Cationic	Ion Exchange	SCX
		Anionic	Ion Exchange	SAX Aminopropyl
		Non-Polar & Cationic	Reversed Phase and Ion Exchange	SOLA SCX/WCX, Verify CX Retain CX
		Non-Polar & Anionic	Reversed Phase and Ion Exchange	SOLA SAX/WAX, Verify AX Retain AX
Aqueous	Organic Soluble	Non-Polar	Reversed Phase	SOLA HRP, C18, C8, Phenyl, Retain PEP
Organic	Organic Soluble	Polar	Normal Phase	Silica, Aminopropyl Cyano, Diol
		Moderately Polar	Normal Phase	Silica, Florisil, Aminopropyl Cyano, Diol
		Cationic	Ion Exchange	SCX
		Anionic	Ion Exchange	SAX Aminopropyl
		Non-Polar & Cationic	Reversed Phase and Ion Exchange	Verify CX
		Non-Polar & Anionic	Reversed Phase and Ion Exchange	Verify AX

Visit our online product selector tool at
www.thermoscientific.com/crc

SPE Phase Selection by Manufacturer

Thermo Scientific™ HyperSep™ Product	Alternative To
HyperSep Retain PEP Page 1-012	Oasis HLB (Waters) Strata-X (Phenomenex) Bond Elut Plexa (Agilent) Easy (Macherey Nagel) ISOLUTE™ ENV+ (Biotage) SampliQ OPT StyreScreen DVB (UCT) H ₂ O-philic DVB (JT Baker)
HyperSep Retain CX Page 1-013	Oasis MCX (Waters) Strata-X-C (Phenomenex) Bond Elut™ Plexa PCX (Agilent) StyreScreen DBX (UCT) H ₂ O-philic SC-DVB (JT Baker)
HyperSep Retain AX Page 1-014	Oasis MAX (Waters) Bond Elut Plexa PAX (Agilent) StyreScreen QAX (UCT) H ₂ O-philic SA-DVB (JT Baker)
HyperSep Hypercarb Page 1-015	SupelClean™ ENVI-Carb™ Bond Elut Carbon
HyperSep C18 Page 1-022	CLEAN-UP™ C18-U Supelclean ENVI-18 / LC-18 BAKERBOND SPE™ Polar Plus™ ISOLUTE C18 CHROMABOND™ C18 Bond Elut C18 Strata C18-U Sep-Pak C18
HyperSep C8 Page 1-023	CLEAN-UP C8 Supelclean ENVI-8 / LC-8 BAKERBOND SPE Octyl C8 ISOLUTE C8 CHROMABOND C8 Bond Elut C8 Strata C8 Sep-Pak C8
HyperSep Phenyl Page 1-024	CLEAN-UP Phenyl Supelclean LC-Ph BAKERBOND SPE Phenyl ISOLUTE Ph Bond Elut Ph Strata Phenyl (PH)
HyperSep Silica Page 1-025	CLEAN-UP Silica Supelclean LC-Si BAKERBOND SPE Silica Gel ISOLUTE Silica CHROMABOND SiOH Bond Elut Si Strata Si-1 Sep-Pak Si

Thermo Scientific HyperSep Product	Alternative To
HyperSep SCX Page 1-026	CLEAN-UP BCX Supelclean LC-SCX BAKERBOND SPE Aromatic Sulfonic Acid ISOLUTE SCX CHROMABOND SA Bond Elut SCX Strata SCX
HyperSep SAX Page 1-027	CLEAN-UP QAX Supelclean LC-SAX BAKERBOND SPE Quaternary amine ISOLUTE SAX CHROMABOND SB Bond Elut SAX Strata SAX
HyperSep Verify CX Page 1-028	CLEAN-UP DAU Discovery DSC-MCAX ISOLUTE HCX CHROMABOND Drug Bond Elut Certify I Strata Screen-C
HyperSep Verify AX Page 1-029	CLEAN-UP THC ISOLUTE HAX Bond Elut Certify II Strata Screen-A
HyperSep Florisil Page 1-030	CLEAN-UP Florisil Supelclean ENVI-Florisil / LC-Florisil BAKERBOND SPE Florisil ISOLUTE Florisil CHROMABOND Florisil Bond Elut Florisil Strata FL-PR Sep-Pak Florisil
HyperSep Aminopropyl Page 1-031	CLEAN-UP Aminopropyl Supelclean LC-NH ₂ BAKERBOND SPE Amino ISOLUTE NH ₂ CHROMABOND NH ₂ Bond Elut NH ₂ Strata NH ₂ Sep-Pak NH ₂
HyperSep Cyano Page 1-032	Supelclean LC-CN BAKERBOND SPE Cyano CHROMABOND CN Bond Elut Cyano Sep-Pak Cyanopropyl Strata CN
HyperSep Diol Page 1-033	Supelclean LC-Diol BAKERBOND SPE Diol CHROMABOND OH Bond Elut Diol Sep-Pak Diol



Removing Uncertainty by Applying the Science to SPE

Our comprehensive range of SPE solutions offer unparalleled performance in purity of extract and reproducibility. Having a fundamental effect on the quality, time and analysis cost, SPE is a critical step during the sample analysis procedure.

We are dedicated to supplying the highest quality SPE solutions, in combination with providing our customers with the support and resources to optimize their SPE solutions and maximize their analysis.

The Importance of SPE

Reducing the effects of the matrix on the separation (GC/LC) and detection (UV, MS etc.) is beneficial. The use of SPE as a sample preparation technique can significantly improve analysis aiding robustness and generating reproducible accurate, precise and sensitive analytical methods. This relies on the ability of SPE to reproducibly:

Maximize detection selectivity

- Reduce ion suppression
- Reduce protein binding
- Reduce matrix interferences by elimination of matrix and particulates
- Compatibility of solvent with analytical technique

Improve analytical system performance

- Longer column lifetimes
- Less detector maintenance
- Autosampler syringes less likely to block
- Less contamination


Improve sensitivity by concentration of the analytes

- Lower limits of detection
- More accurate quantitation
- Improved data processing

Matrix effects are an issue in many areas of analytical chemistry, however, modification of ionization (ion suppression/enhancement) can be a major problem in atmospheric pressure ionization (API) mass spectrometry and, in particular, electrospray based ion sources.



SPE Procedure – Six Steps for a Clean Extract



1. Sample Pre-treatment

It is important to optimize the sample for effective analyte retention. The following should be considered:

- Sample volume/analyte concentration/matrix complexity
- Adjust sample/matrix composition for proper dilution/ionic strength
- Sample pH for optimum retention
- Confirm that analytes are free in solution
- Remove any unwanted particulates via filtration or centrifugation

2. Column Conditioning

Prepare the sorbent for effective interaction(s) with the compounds of interest.

- Use appropriate solvent for column condition/activation
- Prevent sorbent drying during conditioning

3. Column Pre-equilibration

Equilibrate with weakly eluting solvent to prepare the phase for sample addition.

- Use the same solvent as for sample pre-treatment
- Prevent sorbent drying during column equilibration

4. Sample Loading

Analytes are retained on the sorbent.

- Apply samples at appropriate flow rate (1mL/minute typical)

For Reversed-Phase Interactions

- Neutral compounds are not affected by pH
- For charged compounds, a pH at which the compound is not charged is used. Neutralize the molecule according to the following:
 - For basic compounds, the neutral molecule exists at least 2pH units below the pK_a of the compound
 - For acidic compounds, the neutral molecule exists at least 2pH units above the pK_a of the compound

For Normal-Phase Interactions

- pH is not normally an issue in normal phase interactions, as the solvents used are typically non-polar organic solvents, rather than water
- There is no need to verify the sample application pH

For Ion-Exchange Interactions

- pH and pK_a are important considerations
- Acidic compounds are extracted from a sample solution at least 2pH units above the pK_a of the analyte
- Basic compounds are extracted from a sample solution 2 or more pH units below the pK_a of the analyte
- For second (organic) wash, choose the strongest solution where no compound breakthrough occurs
- For elution step, use a solution stronger than where all the compound of interest is eluted
- NB: when choosing these solutions allow some margin for error

5. Wash Away Interferences

Remove impurities bound less strongly than the compounds of interest.

- Select a strong enough wash solvent to remove interferences but weak enough to leave compounds of interest bound
- Selectively rinse away the less strongly bonded interferences
- Wash solvent selected according to phase mechanism/analyte properties

6. Elute Compounds of Interest

Selectively recover the analyte(s) by disrupting the analyte-sorbent interaction.

- Selectively elute analytes of interest using different solvents
- Smaller elution volume produces a more concentrated extract
- Select elution solvent that does not elute strongly retained impurities
- Select elution solvent according to phase mechanism/analyte properties

It is important to optimize the Wash and Elution steps in order to obtain maximum levels of recovery.

Get your copy of our application notebook, featuring over 150 SPE applications
www.thermoscientific.com/spe



Method Development Optimization in SPE


For cleaner extracts, optimization of the SPE process can be important. By optimizing the load, wash, and elution steps of the SPE process, a cleaner sample extract can be obtained.

An example of this is in the development of an SPE method for testosterone.

HyperSep Retain PEP SPE Cartridge


(Part Number **60107-201**)

1. Condition with 1mL methanol followed by 1mL water
2. Load 1mL of 500ng/mL sample in water
3. Sequentially wash with increasing strengths of methanol in water, collecting the eluent 0% methanol/100% water to 90% methanol/10% water, increasing methanol content by 10% each time
4. Elution with multiple volumes of methanol



LogP = 3.6 P_{K_a} = 2.99

Elution Profile of Testosterone



By following this procedure the optimum wash and elution steps can be determined. In subsequent analyses a single wash stage would be used at 50% methanol followed by a single elution step at 90% methanol. This ensures simplicity of the process, optimum recovery and cleanliness of extract.

An alternative would be to use a HyperSep Retain AX SPE cartridge with the added benefit of Phospholipid removal.



For more information on how to optimize your SPE methodology and specific applications visit www.thermoscientific.com/spe and request your free copy of the Thermo Scientific HyperSep Columns Application Notebook.

Thermo Scientific SOLA SPE Cartridges and Plates

Join the revolution

Thermo Scientific™ SOLA™ cartridges and plates provide clean, highly reproducible sample extracts with lower elution volumes due to the unique and innovative frit-less SPE design.

SOLA SPE products provide unparalleled performance characteristics compared to conventional SPE, phospholipid removal and protein precipitation products. This includes:

- Higher levels of reproducibility
- Lower elution volumes
- High extract cleanliness

This increased performance gives higher confidence in analytical results and lowers cost without compromising ease of use or requiring complex method development.

SOLA SPE products provide reduced failure rates, higher analysis speeds and lower solvent requirements which are critical in today's high throughput bio-analytical and clinical research laboratory environments.

Technical Information

The following information highlights the advantages associated with SOLA SPE products over conventional loose packed SPE products.

Improved Reproducibility and Recovery

The data in Figure 1 and Table 1 shows that SOLA SPE products out perform conventional loose packed competitor SPE products even when utilizing the competitor recommended generic methodology.

The information presented shows lower variability sample to sample, and higher recovery levels for SOLA SPE products compared to conventional loose packed competitor SPE products, ensuring you achieve the correct result time after time.




Figure 1: SOLA SPE shows significantly higher reproducibility and recovery levels

	SOLA		Competitor (i)		Competitor (ii)	
	Recovery	Precision (%RSD n=30)	Recovery	Precision (%RSD n=30)	Recovery	Precision (%RSD n=30)
Caffeine	89.88%	4.37%	82.88%	23.89%	63.10%	12.11%
Hydrocortisone	93.76%	3.32%	83.95%	20.49%	76.05%	10.41%
Carbamazepine	98.72%	2.71%	84.47%	21.83%	78.64%	7.79%

Table 1: SOLA SPE provides higher reproducibility and recovery levels

Improved reproducibility

Figure 2 shows that SOLA SPE products have consistent recoveries (across all 30 test samples). The conventional loose packed SPE product from competitor (i) shows that on average one in every four samples gives a significantly lower recovery. This results in inconsistencies in results which potentially would result in costly re-analysis. In comparison SOLA SPE products provide significantly higher levels of reproducibility which is vitally important for high throughput studies.




Figure 2: Shows inconsistency of loose packed products compared to SOLA SPE.



View product information
and application notes

Higher sensitivity and lower solvent consumption

Figure 3 shows that SOLA SPE products achieve excellent recovery levels even with low volumes of extract solvents, resulting in a more concentrated analyte and increased sensitivity. Additional cost and time saving benefits can be achieved from reduced sample dry down time and solvent usage. These low volume extractions would be significantly compromised when using a conventional loose packed, low bed weight, SPE product.






Figure 3: SOLA SPE products recovery and reproducibility at low extraction volumes

Cleanliness of Extract


SOLA SPE products proprietary manufacturing process provides a cleaner product and, as a result, a cleaner sample extract. This is shown in Figure 4 where SOLA SPE products are compared against competitor (i) conventional loose packed SPE product, which have both been extracted with acetonitrile, dichloromethane and methanol, respectively.



Methanol extract comparison:
SOLA products versus competitor (i)



Acetonitrile extract comparison:
SOLA products versus competitor (i)



Dichloromethane extract comparison:
SOLA products versus competitor (i)

Figure 4: Shows SOLA SPE products are significantly cleaner than the equivalent loose packed SPE product from competitor (i)

In today's demanding bio-analytical or clinical research laboratory environment where reproducibility, certainty of results and cost saving are fundamental requirements, SOLA SPE products are an indispensable tool to provide confidence and first time/every time success in the analytical process.



Product Information

SOLA SPE cartridges and plates are available in 10mg/1mL cartridge and 10mg/2mL 96 well plate formats.

SOLA SPE Cartridges

Description	Bed Weight (mg)	Column Volume (mL)	Cat No.	Quantity
SOLA HRP	10	1	60109-001	100 Pack
SOLA SCX	10	1	60109-002	100 Pack
SOLA SAX	10	1	60109-003	100 Pack
SOLA WCX	10	1	60109-004	100 Pack
SOLA WAX	10	1	60109-005	100 Pack

SOLA 96 Well Plates

Description	Bed Weight (mg)	Well Volume (mL)	Cat No.	Quantity
SOLA HRP	10	2	60309-001	1 Pack
SOLA SCX	10	2	60309-002	1 Pack
SOLA SAX	10	2	60309-003	1 Pack
SOLA WCX	10	2	60309-004	1 Pack
SOLA WAX	10	2	60309-005	1 Pack

To learn more about how SOLA SPE can benefit your workflow, method scale down, applications, testimonials and more visit

www.thermoscientific.com/sola-spe



Thermo Scientific HyperSep SPE Cartridges and Plates

HyperSep Retain PEP SPE Cartridges and Plates

Versatile polymeric material for retention of polar and non-polar analytes

- Exceptional recoveries for polar and non-polar analytes
- High and consistent recoveries
- High capacity, high purity, highly porous polystyrene DVB material modified with urea functional groups to give balanced retention of polar and non-polar analytes
- Fast and easy sample preparation and method development
- pH stable 0 to 14



HyperSep Retain PEP SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
30	1	60107-201	100 Pack
30	3	60107-202	50 Pack
60	3	60107-203	50 Pack
60	6	60107-208	30 Pack
100	3	60107-217	50 Pack
100	6	60107-207	30 Pack
150	6	60107-211	30 Pack
200	3	60107-204	50 Pack
200	6	60107-212	30 Pack
500	3	60107-205	50 Pack
500	6	60107-206	30 Pack
1,000	25	60107-215	20 Pack
1,000	6	60107-218	30 Pack
2,000	25	60107-214	20 Pack

Note: For bed weights less than 30mg please refer to SOLA SPE products on page **1-011**.

HyperSep Retain PEP 96 Fixed Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
30	1	60306-207	1 Each
60	1	60306-208	1 Each

Note: For bed weights less than 30mg please refer to SOLA SPE products on page **1-011**.

HyperSep Retain PEP 96 Removable Well Plate

Bed Weight (mg)	Volume (mL)	Cat. No.	Quantity
30	1	60303-207	1 Each

Note: For bed weights less than 30mg please refer to SOLA SPE products on page **1-011**.

Applications:

- Drugs and metabolites in biological fluids
- Peptides in serum, plasma or biological fluids
- Environmental samples

HyperSep Retain CX SPE Cartridges and Plates

Versatile polymeric material for retention of basic compounds

- Exceptional recoveries for basic analytes
- High and consistent recoveries
- High capacity, high purity, highly porous polystyrene DVB material partially modified with sulfonic acid functional groups
- Fast and easy sample preparation and method development
- pH stable 0 to 14

HyperSep Retain CX SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
30	1	60107-301	100 Pack
30	3	60107-302	50 Pack
60	3	60107-303	50 Pack
60	6	60107-308	30 Pack
100	6	60107-307	30 Pack
150	6	60107-311	30 Pack
200	3	60107-304	50 Pack
200	6	60107-314	30 Pack
500	3	60107-305	50 Pack
500	6	60107-306	30 Pack
1,000	25	60107-315	20 Pack
2,000	25	60107-312	20 Pack

Note: For bed weights less than 30mg please refer to SOLA SPE products on page [1-011](#).

HyperSep Retain CX 96 Fixed Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
30	1	60306-303	1 Each
60	1	60306-304	1 Each

Note: For bed weights less than 30mg please refer to SOLA SPE products on page [1-011](#).

HyperSep Retain CX 96 Removable Well Plate

Bed Weight (mg)	Volume (mL)	Cat. No.	Quantity
30	1	60303-307	1 Each

Note: For bed weights less than 30mg please refer to SOLA SPE products on page [1-011](#).

Applications:

- Analysis of a wide range of drugs of abuse including basic and neutral drugs

HyperSep Retain AX SPE Cartridges and Plates

Versatile polymeric material for retention of acidic compounds

- Exceptional recoveries for acidic analytes
- High and consistent recoveries
- High capacity, high purity, highly porous polystyrene DVB material partially modified with quaternary amine functional groups
- Fast and easy sample preparation and method development
- pH stable 0 to 14



HyperSep Retain AX SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
30	1	60107-401	100 Pack
30	3	60107-402	50 Pack
60	3	60107-403	50 Pack
60	6	60107-408	30 Pack
100	3	60107-417	50 Pack
100	6	60107-407	30 Pack
150	6	60107-411	30 Pack
200	3	60107-404	50 Pack
200	6	60107-412	30 Pack
500	3	60107-405	50 Pack
500	6	60107-406	30 Pack
500	15	60107-419	20 Pack
1,000	25	60107-415	20 Pack
2,000	25	60107-414	20 Pack

Note: For bed weights less than 30mg please refer to SOLA SPE products on page [1-011](#).

HyperSep Retain AX 96 Fixed Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
30	1	60306-403	1 Each
60	1	60306-404	1 Each

Note: For bed weights less than 30mg please refer to SOLA SPE products on page [1-011](#).

HyperSep Retain AX 96 Removable Well Plate

Bed Weight (mg)	Volume (mL)	Cat. No.	Quantity
30	1	60303-407	1 Each

Note: For bed weights less than 30mg please refer to SOLA SPE products on page [1-011](#).

Applications:

- Analysis of THC and its metabolites

HyperSep Hypercarb SPE Cartridges and Plates

Unique material for retention of highly polar compounds

- Thermo Scientific™ HyperSep™ Hypercarb™ SPE is a flat, 100% porous graphitic carbon (PGC) with selectivity for structurally similar compounds, offering separation of compounds with simple solvents
- pH stable 0 to 14
- High batch-to-batch reproducibility
- Strong retention properties allow use of low bed weights for concentrated extracts
- Interaction mechanism with polar molecules



HyperSep Hypercarb SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
25	1	60106-304	50 Pack
50	1	60106-303	50 Pack
100	1	60106-302	30 Pack
200	3	60106-301	30 Pack
500	6	60106-402	20 Pack
1,000	6	60106-403	10 Pack
2,000	15	60106-404	10 Pack

HyperSep Hypercarb 96 Fixed Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
10	1	60306-501	1 Each
25	1	60306-502	1 Each
50	1	60306-503	1 Each

HyperSep Servo+ Cartridges and Plates

Thermo Scientific™ HyperSep™ Servo+™ SPE provides optimized SPE solutions for analysis of drugs of abuse. The product range has been specifically designed to provide an out-of-the-box solution for sample preparation of drugs of abuse.


The product range offers high levels of reproducibility and cleanliness of extract for dealing with difficult to analyze biological sample matrices where absolute confidence in your results is required.

Optimized extraction procedures have been developed to allow ready-to-use solutions where there are no requirements for method development and allow confidence that the sample extraction technique will provide the cleanest extract possible.

The following information highlights the elution profiles generated to provide optimized methods for the different drug classes for Servo+ SPE.

HyperSep Servo+ Opiates SPE

Elution Profile




Optimized Method:

Washing the Servo+ Opiates Cartridge or plate with 2% Formic Acid in Methanol, followed by 10:90 Methanol:Water/ NH_4OH (95:5) provides an extract with minimal interferences.

HyperSep Servo+ AMP SPE

Elution Profile




Optimized Method:

Washing the Servo+ AMP cartridge with 2% Formic Acid in Methanol, followed by 80:20 Methanol:Water/ NH_4OH (95:5) provides an extract with minimal interferences.

HyperSep Servo+ PCP SPE

Elution Profile




Optimized Method:

Washing the Servo+ PCP cartridge or plate with 2% Formic Acid in Methanol, followed by 80:20 Methanol:Water/ NH_4OH (95:5) provides an extract with minimal interferences.

HyperSep Servo+ Cocaine SPE

Elution Profile




Optimized Method:

Washing the Servo+ Cocaine cartridge or plate with 2% Formic Acid in Methanol, followed by 50:50 Methanol:Water/ NH_4OH (95:5) provides an extract with minimal interferences..

HyperSep Servo+ THC SPE

Elution Profile



Optimized Method:

Washing the Servo+ THC cartridge or plate with 2% Ammonia hydroxide in Methanol, followed by 70:30 Methanol:Water/Formic Acid (95:5) provides an extract with minimal interferences.

HyperSep Servo+ Total SPE Cartridges

Servo+ Total – for screening of a wide range of unknown components

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
30	1	60110-101C	100 Pack
30	3	60110-102C	50 Pack
60	3	60110-103C	50 Pack
200	3	60110-104C	50 Pack
500	3	60110-105C	50 Pack
500	6	60110-106C	30 Pack
100	6	60110-107C	30 Pack
60	6	60110-108C	30 Pack
150	6	60110-109C	30 Pack
200	6	60110-110C	30 Pack
1,000	25	60110-111C	20 Pack

Note: For bed weights less than 30mg please refer to SOLA SPE products on page [1-011](#).

HyperSep Servo+ Total 96 Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
30	1	60110-103P	1 Each

Note: For bed weights less than 30mg please refer to SOLA SPE products on page [1-011](#).

HyperSep Servo+ Total-A SPE Cartridges

Servo+ Total A – for screening of a wide range of unknown acidic compounds

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
30	1	60110-201C	100 Pack
30	3	60110-202C	50 Pack
60	3	60110-203C	50 Pack
200	3	60110-204C	50 Pack
500	3	60110-205C	50 Pack
500	6	60110-206C	30 Pack
100	6	60110-207C	30 Pack
60	6	60110-208C	30 Pack
150	6	60110-209C	30 Pack
200	6	60110-210C	30 Pack
1,000	25	60110-211C	20 Pack

Note: For bed weights less than 30mg please refer to SOLA SPE products on page [1-011](#).

HyperSep Servo+ Total-A 96 Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
30	1	60110-203P	1 Each

Note: For bed weights less than 30mg please refer to SOLA SPE products on page [1-011](#).

HyperSep Servo+ Total-B SPE Cartridges

Servo+ Total B – for screening of a wide range of unknown basic compounds

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
30	1	60110-301C	100 Pack
30	3	60110-302C	50 Pack
60	3	60110-303C	50 Pack
200	3	60110-304C	50 Pack
500	3	60110-305C	50 Pack
500	6	60110-306C	30 Pack
100	6	60110-307C	30 Pack
60	6	60110-308C	30 Pack
150	6	60110-309C	30 Pack
200	6	60110-310C	30 Pack
1,000	25	60110-311C	20 Pack

Note: For bed weights less than 30mg please refer to SOLA SPE products on page **1-011**.

HyperSep Servo+ Total-B 96 Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
30	1	60110-303P	1 Each

Note: For bed weights less than 30mg please refer to SOLA SPE products on page **1-011**.

HyperSep Servo+ THC SPE Cartridges

Servo+ THC – for dedicated extraction of THC and associated metabolites

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
30	1	60110-401C	100 Pack
30	3	60110-402C	50 Pack
60	3	60110-403C	50 Pack
200	3	60110-404C	50 Pack
500	3	60110-405C	50 Pack
500	6	60110-406C	30 Pack
100	6	60110-407C	30 Pack
60	6	60110-408C	30 Pack
150	6	60110-409C	30 Pack
200	6	60110-410C	30 Pack
1,000	25	60110-411C	20 Pack

Note: For bed weights less than 30mg please refer to SOLA SPE products on page **1-011**.

HyperSep Servo+ THC 96 Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
30	1	60110-403P	1 Each

Note: For bed weights less than 30mg please refer to SOLA SPE products on page **1-011**.

HyperSep Servo+ OPIATE SPE Cartridges

Servo+ Opiates – for dedicated extraction of opiate class compounds

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
30	1	60110-501C	100 Pack
30	3	60110-502C	50 Pack
60	3	60110-503C	50 Pack
200	3	60110-504C	50 Pack
500	3	60110-505C	50 Pack
500	6	60110-506C	30 Pack
100	6	60110-507C	30 Pack
60	6	60110-508C	30 Pack
150	6	60110-509C	30 Pack
200	6	60110-510C	30 Pack
1,000	25	60110-511C	20 Pack

Note: For bed weights less than 30mg please refer to SOLA SPE products on page [1-011](#).

HyperSep Servo+ OPIATE 96 Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
30	1	60110-503P	1 Each

Note: For bed weights less than 30mg please refer to SOLA SPE products on page [1-011](#).

HyperSep Servo+ PCP SPE Cartridges

Servo+ PCP – for dedicated extraction of PCP class compounds

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
30	1	60110-601C	100 Pack
30	3	60110-602C	50 Pack
60	3	60110-603C	50 Pack
200	3	60110-604C	50 Pack
500	3	60110-605C	50 Pack
500	6	60110-606C	30 Pack
100	6	60110-607C	30 Pack
60	6	60110-608C	30 Pack
150	6	60110-609C	30 Pack
200	6	60110-610C	30 Pack
1,000	25	60110-611C	20 Pack

Note: For bed weights less than 30mg please refer to SOLA SPE products on page [1-011](#).

HyperSep Servo+ PCP 96 Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
30	1	60110-603P	1 Each

Note: For bed weights less than 30mg please refer to SOLA SPE products on page [1-011](#).

HyperSep Servo+ AMP SPE Cartridges

Servo+ AMP – for dedicated extraction of amphetamine class compounds

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
30	1	60110-701C	100 Pack
30	3	60110-702C	50 Pack
60	3	60110-703C	50 Pack
200	3	60110-704C	50 Pack
500	3	60110-705C	50 Pack
500	6	60110-706C	30 Pack
100	6	60110-707C	30 Pack
60	6	60110-708C	30 Pack
150	6	60110-709C	30 Pack
200	6	60110-710C	30 Pack
1,000	25	60110-711C	20 Pack

Note: For bed weights less than 30mg please refer to SOLA SPE products on page **1-011**.

HyperSep Servo+ AMP 96 Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
30	1	60110-703P	1 Each

Note: For bed weights less than 30mg please refer to SOLA SPE products on page **1-011**.

HyperSep Servo+ Cocaine SPE Cartridges

Servo+ Cocaine – for dedicated extraction of cocaine class compounds

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
30	1	60110-801C	100 Pack
30	3	60110-802C	50 Pack
60	3	60110-803C	50 Pack
200	3	60110-804C	50 Pack
500	3	60110-805C	50 Pack
500	6	60110-806C	30 Pack
100	6	60110-807C	30 Pack
60	6	60110-808C	30 Pack
150	6	60110-809C	30 Pack
200	6	60110-810C	30 Pack
1,000	25	60110-811C	20 Pack

Note: For bed weights less than 30mg please refer to SOLA SPE products on page **1-011**.

HyperSep Servo+ Cocaine 96 Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
30	1	60110-803P	1 Each

Note: For bed weights less than 30mg please refer to SOLA SPE products on page **1-011**.

To learn more about HyperSep Servo+ visit
www.thermoscientific.com/servo



HyperSep C18 SPE Cartridges and Plates

Feature a highly retentive alkyl-bonded phase for non-polar to moderately polar compounds

- Retentive for most non-polar compounds
- Retains most organic analytes from aqueous matrices
- Hydrophobic reversed phase
- For non-polar to moderately polar compounds

HyperSep C18 SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
25	1	60108-376	100 Pack
50	1	60108-390	100 Pack
100	1	60108-302	100 Pack
100	3	60108-765	50 Pack
200	3	60108-303	50 Pack
500	3	60108-304	50 Pack
500	6	60108-305	30 Pack
500	10	60108-786	50 Pack
1,000	6	60108-301	30 Pack
2,000	15	60108-701	20 Pack
5,000	25	60108-702	20 Pack
10,000	75	60108-703	10 Pack

Applications:

- Retentive for non-polar compounds
- Retains most organic analytes from aqueous matrices

HyperSep C18 96 Fixed Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
10	1	60307-201	1 Each
25	1	60307-202	1 Each
50	1	60307-203	1 Each
100	1	60307-204	1 Each

HyperSep C18 96 Removable Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
10	1	60300-425	1 Each
25	1	60300-426	1 Each
50	1	60300-427	1 Each
100	1	60300-428	1 Each

HyperSep C8 SPE Cartridges and Plates

Less retentive alternative to C18 for polar and non-polar compounds

- Well-suited for methods requiring less retention than C18
- Excellent for moderately polar analytes
- Potential for polar interactions greater than with C18



HyperSep C8 SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
50	1	60108-391	100 Pack
100	1	60108-392	100 Pack
200	3	60108-393	50 Pack
500	3	60108-309	50 Pack
500	6	60108-394	30 Pack
1,000	6	60108-427	30 Pack
2,000	15	60108-704	20 Pack
5,000	25	60108-705	20 Pack
10,000	75	60108-706	10 Pack

HyperSep C8 96 Fixed Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
10	1	60307-211	1 Each
25	1	60307-212	1 Each
50	1	60307-213	1 Each
100	1	60307-214	1 Each

HyperSep C8 96 Removable Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
10	1	60300-445	1 Each
25	1	60300-446	1 Each
50	1	60300-447	1 Each
100	1	60300-448	1 Each

Applications:

- Drugs and their metabolites in biological samples
- Peptides in biological samples

HyperSep Phenyl SPE Cartridges and Plates

Offer alternative selectivity for retention of basic compounds

- Commonly employed for non-polar interactions
- The electron density of the aromatic ring provides different selectivity than other non-polar sorbents

HyperSep Phenyl SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
50	1	60108-516	100 Pack
100	1	60108-386	100 Pack
200	3	60108-387	50 Pack
500	3	60108-388	50 Pack
500	6	60108-389	30 Pack
1,000	6	60108-517	30 Pack
2,000	15	60108-707	20 Pack
5,000	25	60108-708	20 Pack
10,000	75	60108-709	10 Pack

Applications:

- Extraction of aromatic compounds
- Extraction of basic analytes

HyperSep Phenyl 96 Fixed Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
10	1	60307-221	1 Each
25	1	60307-222	1 Each
50	1	60307-223	1 Each
100	1	60307-224	1 Each

HyperSep Phenyl 96 Removable Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
10	1	60300-685	1 Each
25	1	60300-686	1 Each
50	1	60300-687	1 Each
100	1	60300-688	1 Each

HyperSep Silica SPE Cartridges and Plates

A polar sorbent primarily used to retain analytes in non-polar matrices

- Effectively separates compounds of very similar structure
- Adsorbs analytes from non-polar solvents like hydrocarbons, less polar esters and ethers
- Suitable for use as an intermediate strength cation exchanger in aqueous medium

HyperSep Silica SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
50	1	60108-409	100 Pack
100	1	60108-317	100 Pack
200	3	60108-410	50 Pack
500	3	60108-315	50 Pack
500	6	60108-411	30 Pack
500	10	60108-793	50 Pack
1,000	6	60108-426	30 Pack
2,000	15	60108-710	20 Pack
5,000	25	60108-711	20 Pack
10,000	25	60108-853	20 Pack
10,000	75	60108-712	10 Pack
50,000	150	60108-850	10 Pack
20,000	75	60108-851	10 Pack
70,000	150	60108-852	10 Pack

HyperSep Silica 96 Fixed Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
10	1	60307-231	1 Each
25	1	60307-232	1 Each
50	1	60307-233	1 Each
100	1	60307-234	1 Each

HyperSep Silica 96 Removable Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
10	1	60300-485	1 Each
25	1	60300-486	1 Each
50	1	60300-487	1 Each
100	1	60300-488	1 Each

Applications:

- Extraction of polar compounds including aldehydes, amines, drugs, pesticides and herbicides
- Extraction of carotenoids, fat-soluble vitamins, aflatoxins in food matrices
- Extraction of fatty acids and phospholipids

HyperSep SCX Strong Cation Exchanger SPE Cartridges and Plates

A strong cation exchange sorbent for extraction of charged basic compounds

- Modified benzene ring provides potential for non-polar interactions
- Low pK_a value sorbent

HyperSep SCX SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
50	1	60108-420	100 Pack
100	1	60108-421	100 Pack
200	3	60108-422	50 Pack
500	3	60108-423	50 Pack
500	6	60108-520	30 Pack
1,000	6	60108-433	30 Pack
2,000	15	60108-716	20 Pack
5,000	25	60108-717	20 Pack
10,000	75	60108-718	10 Pack
20,000	75	60108-857	10 Pack
50,000	150	60108-858	10 Pack

HyperSep SCX 96 Fixed Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
10	1	60307-251	1 Each
25	1	60307-252	1 Each
50	1	60307-253	1 Each
100	1	60307-254	1 Each

HyperSep SCX 96 Removable Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
10	1	60300-585	1 Each
25	1	60300-586	1 Each
50	1	60300-587	1 Each
100	1	60300-588	1 Each

Applications:

- Isolation of anionic proteins
- Removal of acidic food pigments
- Isolation of phenolic compounds
- Nucleic acids, nucleotides and surfactants

HyperSep SAX Strong Anion Exchanger SPE Cartridges and Plates

Strong anion exchange sorbent for extraction of weak acids

- Extracts negatively charged compounds from both aqueous and non-aqueous solutions
- Ideally suited for the extraction of weak acids, e.g. carboxylic acids
- Amine group masks the effect of the carbon chain of the functional group
- Selectivity can be tuned by selection of buffer

HyperSep SAX SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
50	1	60108-417	100 Pack
100	1	60108-418	100 Pack
200	3	60108-419	50 Pack
500	3	60108-521	50 Pack
500	6	60108-360	30 Pack
1,000	6	60108-434	30 Pack
2,000	15	60108-713	20 Pack
5,000	25	60108-714	20 Pack
10,000	75	60108-715	10 Pack

HyperSep SAX 96 Fixed Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
10	1	60307-241	1 Each
25	1	60307-242	1 Each
50	1	60307-243	1 Each
100	1	60307-244	1 Each

HyperSep SAX 96 Removable Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
10	1	60300-565	1 Each
25	1	60300-566	1 Each
50	1	60300-567	1 Each
100	1	60300-568	1 Each

Applications:

- Isolation of anionic proteins
- Removal of acidic food pigments
- Isolation of phenolic compounds
- Nucleic acids, nucleotides and surfactants

HyperSep Verify CX SPE Cartridges and Plates

Features non-polar and ionic characteristics for improved analysis of drugs of abuse

- Thermo Scientific™ HyperSep™ Verify™ CX mixed mode sorbent for extraction of basic drugs
- Separation based on two functional groups: reversed phase C8 and an ion exchanger, benzene sulfonic acid
- Co-polymerized on a rigid, purified silica gel support

HyperSep Verify CX SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
50	1	60108-741	100 Pack
130	1	60108-719	100 Pack
130	10	60108-769	50 Pack
200	3	60108-777	50 Pack
200	6	60108-722	50 Pack
200	10	60108-742	50 Pack
300	3	60108-720	50 Pack
500	3	60108-721	50 Pack
500	6	60108-723	30 Pack
1,000	6	60108-724	30 Pack

HyperSep Verify CX 96 Fixed Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
10	1	60307-261	1 Each
25	1	60307-262	1 Each
50	1	60307-263	1 Each
100	1	60307-264	1 Each

HyperSep Verify CX 96 Removable Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
10	1	60300-805	1 Each
25	1	60300-806	1 Each
50	1	60300-807	1 Each
100	1	60300-808	1 Each

Applications:

- Analysis of a wide range of drugs of abuse from biological matrices, including basic and neutral drugs

HyperSep Verify AX SPE Cartridges and Plates

Features non-polar and ionic characteristics for improved analysis of acidic drugs and metabolites

- Mixed mode sorbent suited for the extraction of acidic drugs and metabolites from biological matrices
- Separation based on two functional groups: reversed phase C8 and an ion exchanger, quaternary amine
- Rigid, purified silica gel support to which the two functionalities are co-polymerized



HyperSep Verify AX SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
130	1	60108-727	100 Pack
130	10	60108-767	50 Pack
200	6	60108-730	50 Pack
200	10	60108-764	50 Pack
300	3	60108-728	50 Pack
500	3	60108-729	50 Pack
500	6	60108-731	30 Pack
1,000	6	60108-732	30 Pack

HyperSep Verify AX 96 Fixed Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
10	1	60307-271	1 Each
25	1	60307-272	1 Each
50	1	60307-273	1 Each
100	1	60307-274	1 Each

HyperSep Verify AX 96 Removable Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
10	1	60300-813	1 Each
25	1	60300-814	1 Each
50	1	60300-815	1 Each
100	1	60300-816	1 Each

Applications:

- Analysis of THC and its metabolites

HyperSep Florisil SPE Cartridges and Plates

Ideal for the isolation of polar compounds from non-polar matrices

- Magnesia-loaded silica gel offers rapid flow rates for large sample volumes
- Extremely polar in nature, especially well-suited for extraction of chlorinated pesticides



HyperSep Florisil SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
50	1	60108-402	100 Pack
100	1	60108-403	100 Pack
200	3	60108-404	50 Pack
500	3	60108-405	50 Pack
500	6	60108-500	30 Pack
1,000	6	60108-431	30 Pack
2,000	15	60108-735	20 Pack
5,000	25	60108-736	20 Pack
10,000	75	60108-737	10 Pack

HyperSep Florisil 96 Fixed Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
10	1	60307-281	1 Each
25	1	60307-282	1 Each
50	1	60307-283	1 Each
100	1	60307-284	1 Each

HyperSep Florisil 96 Removable Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
10	1	60300-725	1 Each
25	1	60300-726	1 Each
50	1	60300-727	1 Each
100	1	60300-728	1 Each

Applications:

- Extraction of pesticides using official methods
- Polychlorinated biphenyls in transformer oil
- Alcohol, aldehydes, amines and drugs

HyperSep Aminopropyl SPE Cartridges and Plates

A polar sorbent for both polar and anion exchange interactions

- Weak anion exchanger
- Enhanced retention of strong anions like sulfonic acids
- Excellent retention of drugs, metabolites and structural isomer

HyperSep Aminopropyl SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
50	1	60108-424	100 Pack
100	1	60108-364	100 Pack
200	3	60108-425	50 Pack
500	3	60108-518	50 Pack
500	6	60108-519	30 Pack
1,000	6	60108-432	30 Pack
2,000	15	60108-738	20 Pack
5,000	25	60108-739	20 Pack
10,000	75	60108-740	10 Pack

HyperSep Aminopropyl 96 Fixed Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
10	1	60307-291	1 Each
25	1	60307-292	1 Each
50	1	60307-293	1 Each
100	1	60307-294	1 Each

HyperSep Aminopropyl 96 Removable Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
10	1	60300-505	1 Each
25	1	60300-506	1 Each
50	1	60300-507	1 Each
100	1	60300-508	1 Each

Applications:

- Separation of structural isomers
- Drugs and metabolites in biological fluids
- Separation of saccharides, phenols and petroleum products

HyperSep Cyano SPE Cartridges and Plates

Optimized for the retention of polar compounds from non-polar matrices

- Medium polarity sorbent, less retentive than silica or diol
- Suitable for reverse phase extractions of moderately polar compounds

HyperSep Cyano SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
50	1	60108-746	100 Pack
100	1	60108-745	100 Pack
200	3	60108-747	50 Pack
500	3	60108-748	50 Pack
500	6	60108-749	30 Pack
1,000	6	60108-750	30 Pack
2,000	15	60108-751	20 Pack
5,000	25	60108-752	20 Pack
10,000	75	60108-753	10 Pack

HyperSep Cyano 96 Fixed Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
10	1	60307-301	1 Each
25	1	60307-302	1 Each
50	1	60307-303	1 Each
100	1	60307-304	1 Each

HyperSep Cyano 96 Removable Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
10	1	60300-821	1 Each
25	1	60300-822	1 Each
50	1	60300-823	1 Each
100	1	60300-824	1 Each

Applications:

- Retaining polar compounds from hexane and oils
- Reversed phase extraction of moderately polar compounds

HyperSep Diol SPE Cartridges and Plates

Ideal for the extraction of polar compounds

- Short alkyl chains with polar functional groups

HyperSep Diol SPE Cartridges

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
50	1	60108-571	100 Pack
100	1	60108-572	100 Pack
200	3	60108-573	50 Pack
500	3	60108-574	50 Pack
500	6	60108-575	30 Pack
1,000	6	60108-576	30 Pack
2,000	15	60108-755	20 Pack
5,000	25	60108-756	20 Pack
10,000	75	60108-757	10 Pack

HyperSep Diol 96 Fixed Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
10	1	60307-311	1 Each
25	1	60307-312	1 Each
50	1	60307-313	1 Each

HyperSep Diol 96 Removable Well Plates

Bed Weight (mg)	Well Volume (mL)	Cat. No.	Quantity
10	1	60300-630	1 Each
25	1	60300-631	1 Each
50	1	60300-632	1 Each
100	1	60300-633	1 Each

Applications:

- Normal phase extraction of polar compounds
- Purification of polar compounds

HyperSep SPE Accessories

A range of accessories to complement the HyperSep SPE product range

HyperSep SPE Accessories

Description	Cat. No.	Quantity
Base Plate for HyperSep-96 Well Plate	60300-301	1 Each
Base Plate for HyperSep-96 Well Plate	60300-303	5 Pack
Sample Collection Plate, 1mL	60300-402	50 Pack
Sample Collection Plate, 2mL	60300-403	50 Pack
Adaptors for 1mL, 3mL and 6mL SPE Cartridges	60104-259	15 Pack
Empty 1mL Wells	60300-318	100 Pack
Empty 1mL Wells, Fritted	60300-311	100 Pack



HyperSep Protein Precipitation Plate

Provides a quick effective approach for removal of proteins from biological compounds using the protein crash technique. In combination with SPE and SLE, the Protein Precipitation Plate offers a comprehensive range of options for sample preparation of biological based compounds. The 96 well plate format provides a high throughput functionality which lends itself to full automation.

- Dual frit design
- Hydrophobic/oleophobic frits to enable only precipitation of proteins
- Pore size optimized for ideal flow rate
- Specially selected polypropylene for low extractables



HyperSep Protein Precipitation Plate

Description	Cat. No.	Quantity
Protein Precipitation Plate	60304-201	1 Each

HyperSep SLE Cartridges and Plates

Solid supported liquid/liquid extraction (SLE) is a fast effective sample preparation technique that provides considerable benefits over liquid-liquid extraction (LLE) protocols for removal of phospholipids from biological samples.

SLE offers the following advantages:

- Greater reproducibility and recoveries compared to LLE techniques
- Prevents emulsification often associated with LLE as the sample and water immiscible solvents are not in direct contact
- Reduced solvent requirements compared to LLE
- Can be completely automated unlike LLE
- Improved cleanliness of sample extract compared to protein precipitation techniques
- Improved sensitivity compared to protein precipitation techniques

HyperSep SLE is available in cartridge and 96 well plate formats and in acidified or basified formats to help simplify the protocol when acidic or basic mobile phases are used.



HyperSep SLE Cartridges

Special treated diatomite SLE (pH=7)

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
200	3	60109-200-3-7	50 Pack
500	3	60109-500-3-7	50 Pack
500	6	60109-500-6-7	30 Pack
1,000	6	60109-1000-6-7	30 Pack
2,000	12	60109-2000-12-7	20 Pack
4,000	25	60109-4000-25-7	15 Pack
20,000	60	60109-20000-60-7	10 Pack

HyperSep SLE 96 Well Plates

Special treated diatomite SLE (pH=7)

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
200	2	60109-200-2-7W	1 Each
300	2	60109-300-2-7W	1 Each
400	2	60109-400-2-7W	1 Each
500	2	60109-500-2-7W	1 Each

HyperSep SLE Cartridges

Special treated diatomite SLE (pH=9)

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
200	3	60109-200-3-9	50 Pack
500	3	60109-500-3-9	50 Pack
500	6	60109-500-6-9	30 Pack
1,000	6	60109-1000-6-9	30 Pack
2,000	12	60109-2000-12-9	20 Pack
4,000	25	60109-4000-25-9	15 Pack
20,000	60	60109-20000-60-9	10 Pack

HyperSep SLE 96 Well Plates

Special treated diatomite SLE (pH=9)

Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
200	2	60109-200-2-9W	1 Each
300	2	60109-300-2-9W	1 Each
400	2	60109-400-2-9W	1 Each
500	2	60109-500-2-9W	1 Each

To learn more visit www.thermoscientific.com/sle

HyperSep Lab Plates

For the purification and sample preparation of proteins, DNA, RNA and other biomolecules

- Sample concentration of small-scale samples
- Available in a range of chromatographic materials
- 96 well plate format with media embedded at the bottom of the plate
- Can be processed manually or using a liquid-handling robot
- Not suitable for use with a vacuum

HyperSep Lab Plates

Description	Cat. No.	Quantity
Polystyrene		
C18	60110-201	5 Pack
C8	60110-202	5 Pack
C4	60110-203	5 Pack
Hypercarb	60110-204	5 Pack
Hypercarb and C18 (Mixed Mode)	60110-205	5 Pack
Zirconium Dioxide	60110-206	5 Pack
Titanium Dioxide	60110-207	5 Pack
SCX	60110-208	5 Pack
SAX	60110-209	5 Pack
Polypropylene		
C18	60110-301	5 Pack
C8	60110-302	5 Pack
C4	60110-303	5 Pack
Hypercarb	60110-304	5 Pack
Hypercarb and C18 (Mixed Mode)	60110-305	5 Pack
Zirconium Dioxide	60110-306	5 Pack
Titanium Dioxide	60110-307	5 Pack
SCX	60110-308	5 Pack
SAX	60110-309	5 Pack

Applications:

- Tissue culture and separation of products
- Sample concentration
- Sample clean-up
- Collection of sample after chromatography

HyperSep Filter Plates

For effective clean-up of small-scale samples

- 96 Well Plate for the purification and separation of proteins, peptides, DNA, RNA and other biomolecules
- Cleanup of microgram-level samples
- Can be used under vacuum
- Available in a range of chromatographic materials

HyperSep Filter Plates

Description	Cat. No.	Quantity
5-7µL Bed Volume		
C18	60110-401	1 Each
C8	60110-402	1 Each
C4	60110-403	1 Each
Hypercarb	60110-404	1 Each
Hypercarb and C18 (Mix Mode)	60110-405	1 Each
Zirconium Dioxide	60110-406	1 Each
Titanium Dioxide	60110-407	1 Each
SCX	60110-408	1 Each
SAX	60110-409	1 Each
40µL Bed Volume		
C18	60110-501	1 Each
C8	60110-502	1 Each
C4	60110-503	1 Each
Hypercarb	60110-504	1 Each
Hypercarb and C18 (Mix Mode)	60110-505	1 Each
Zirconium Dioxide	60110-506	1 Each
Titanium Dioxide	60110-507	1 Each
SCX	60110-508	1 Each
SAX	60110-509	1 Each

Dionex OnGuard II Cartridges

- Thermo Scientific™ Dionex™ OnGuard II™ cartridges remove matrix interferences such as phenolics compounds, metals, cations, anions, or hydrophobic substances encountered in many IC and HPLC applications
- Cartridges are available in 1mL and 2.5mL high-capacity formats



Related Literature

For detailed specifications and applications, see the following PDF documents under Literature on www.thermoscientific.com/dionex

Product Data Sheets

Dionex OnGuard II Sample Pretreatment Cartridges and Workstation

Application Notes

- AN 101: Trace Level Determination of Bromate in Ozonated Drinking Water Using Ion Chromatography
- AN 136: Determination of Inorganic Oxyhalide Disinfection Byproduct Anions and Bromide in Drinking Water Using Ion Chromatography with the Addition of a Postcolumn Reagent for Trace Bromate Analysis
- AN 37: The Determination of Iodide in Milk Products

Application Updates

- AU 140: The Determination of Iodide in Urine

Dionex OnGuard II A Cartridges

- Dionex OnGuard II A cartridges are used in the removal of anionic contaminants from sample matrices and for the neutralization of highly acidic samples
- Cartridges contain styrene-based, anion exchange resin in the bicarbonate form

Key Specifications

- Functionality: anion-exchange, bicarbonate form
- Capacity ($\mu\text{eq}/\text{cartridge}$): 0.7 (1mL cartridge); 1.75 (2.5mL cartridge)
- Solvents: 0-100% HPLC
- pH stability: 0-14
- Mode: removal of anions; pH adjustments of acidic samples

Dionex OnGuard II A Cartridges

Description	Volume (mL)	Cat. No.	Quantity
Dionex OnGuard II A Cartridges	1	057091	48 Pack
Dionex OnGuard II A Cartridges	2.5	057092	48 Pack

Dionex OnGuard II Ag Cartridges

- The Dionex OnGuard II Ag cartridges contain a silver-form, high-capacity, sulfonated, cation-exchange resin similar to the OnGuard II H packing
- Cartridges remove chloride, bromide, and iodide from sample matrices
- Dionex OnGuard II H cartridges should be used after Dionex OnGuard II Ag cartridges to remove dissolved Ag⁺

Key Specifications

- Functionality: cation-exchange, silver form
- Capacity (µeq/cartridge): 2.0-2.2 (1mL); 5.0-5.5 (2.5mL)
- Solvents: 0-100% HPLC
- pH Stability: 0-14
- Mode: removal of chloride, bromide, iodide by precipitation

Dionex OnGuard II Ag Cartridges

Description	Volume (mL)	Cat. No.	Quantity
Dionex OnGuard II Ag Cartridges	1	057089	48 Pack
Dionex OnGuard II Ag Cartridges	1	082762	12 Pack
Dionex OnGuard II Ag Cartridges	2.5	057090	48 Pack

Dionex OnGuard II Ag/H Cartridges

The Dionex OnGuard II Ag/H is a layered cartridge that contains both OnGuard II Ag and Dionex OnGuard II H resins.

- Easily removes chloride, bromide, and iodide from concentrated matrices such as brines
- Traps soluble silver and other cations
- Removes high levels of alkaline earth and transition metals
- Neutralizes caustic samples
- Removes carbonate
- The two-layer cartridge replaces two cartridges in series, and provides greater silver capacity

Dionex OnGuard II Ag/H Cartridges

Description	Volume (mL)	Cat. No.	Quantity
Dionex OnGuard II Ag-H Cartridges	2.5	057410	48 Pack
Dionex OnGuard II Ag-H Cartridges	2.5	082756	12 Pack

Dionex OnGuard II Ba Cartridges

- The Dionex OnGuard II Ba resin is a styrene-based, sulfonic acid resin in the barium form, designed for the removal of high concentrations of sulfate from sample matrices
- For reproducible, quantitative determinations in low-ionic strength samples, activate these cartridges by adding sodium chloride or other sodium salt
- Samples treated with NaCl should be passed through a Dionex OnGuard II Ag cartridge to remove the added chloride, followed by the Dionex OnGuard II H cartridge or Thermo Scientific™ MetPac™ CC-1 chelating column to remove residual silver counterions

Key Specifications

- Functionality: cation-exchange, barium form
- Capacity (µeq/cartridge): 2.0-2.2 (1mL); 5.0-5.5 (2.5mL)
- Solvents: 0-100% HPLC
- pH Stability: 0-14
- Mode: removal of sulfate by precipitation

Dionex OnGuard II Ba Cartridges

Description	Volume (mL)	Cat. No.	Quantity
Dionex OnGuard II Ba Cartridges	1	057093	48 Pack
Dionex OnGuard II Ba Cartridges	1	082763	12 Pack
Dionex OnGuard II Ba Cartridges	2.5	057094	48 Pack

Dionex OnGuard II Ba/Ag/H Cartridges

The Dionex OnGuard II Ba/Ag/H is a layered cartridge containing Dionex OnGuard II Ba, Ag, and H styrene-based, sulfonic acid resins.

- The Ba resin removes high concentrations of sulfate from sample matrices
- The Ag form easily removes chloride, bromide, and iodide from concentrated matrices
- The H form is highly selective for polyvalent cations such as calcium and transition metals
- Ideal for the removal of high levels of alkaline earth and transition metals from sample matrices, neutralization of caustic samples, and removal of carbonate
- Three-layer cartridge can be used in place of three single cartridges in series and has the added advantage of higher silver capacity

Dionex OnGuard II Ba/Ag/H Cartridges

Description	Volume (mL)	Cat. No.	Quantity
Dionex OnGuard II Ba/Ag/H Cartridges	2.5	063955	48 Pack
Dionex OnGuard II Ba/Ag/H Cartridges	2.5	082764	12 Pack

Dionex OnGuard II H Cartridges

- The Dionex OnGuard II H cartridge removes high levels of alkaline earth and transition metals from sample matrices and neutralizes highly alkaline samples such as sodium hydroxide or sodium carbonate. Carbonate can then be removed by sparging the sample.
- Cartridges contain 16% crosslinked, styrene-based, sulfonic acid resin in the hydrogen form. This resin is designed to have very high selectivity for multivalent cations such as calcium and transition metals.

Key Specifications

- Functionality: cation-exchange hydronium form
- Capacity ($\mu\text{eq}/\text{cartridge}$): 2.0-2.2 (1mL); 5.0-5.5 (2.5mL)
- Solvents: 0-100% HPLC
- pH Stability: 0-14
- Mode: removal of alkaline earth and transition metals; pH adjustment of basic samples

Dionex OnGuard II H Cartridges

Description	Volume (mL)	Cat. No.	Quantity
Dionex OnGuard II H Cartridges	1	057085	48 Pack
Dionex OnGuard II H Cartridges	1	082761	12 Pack
Dionex OnGuard II H Cartridges	2.5	057086	48 Pack

Dionex OnGuard II M Cartridges

- The Dionex OnGuard II M cartridge is used for the removal of transition metals and for matrix elimination of alkali and alkaline earth metals
- Cartridges contain an iminodiacetate resin in the ammonium form, ready to use with no lengthy preparation required

Key Specifications

- Functionality: iminodiacetate, ammonium form
- Capacity ($\mu\text{eq}/\text{cartridge}$): 0.4 (1mL); 1.0 (2.5mL)
- Solvents: 0-100% HPLC
- pH Stability: 0-14 (resin shrinks in acid form)
- Mode: concentration of transition metals by chelation (2.5mL format); removal of transition metals (1mL format)

Dionex OnGuard II M Cartridges

Description	Volume (mL)	Cat. No.	Quantity
Dionex OnGuard II M Cartridges	1	057137	48 Pack
Dionex OnGuard II M Cartridges	2.5	057095	48 Pack

Dionex OnGuard II Na Cartridges

- The Dionex OnGuard II Na cartridge removes high levels of alkaline earth and transition metals from sample matrices without acidifying the sample. This ensures good recovery of acid-labile analytes such as nitrite.
- Cartridges contain 16% crosslinked, styrene-based, sulfonic acid resin in the sodium form
- This resin has very high selectivity for multivalent cations such as calcium, magnesium and transition metals

Key Specifications

- Functionality: cation-exchange, sodium form
- Capacity ($\mu\text{eq}/\text{cartridge}$): 2.0-2.2 (1mL); 5.0-5.5 (2.5mL)
- Solvents: 0-100% HPLC
- pH Stability: 0-14
- Mode: removal of alkaline earth and transition metals without acidifying the sample

Dionex OnGuard II Na Cartridges

Description	Volume (mL)	Cat. No.	Quantity
Dionex OnGuard II Na Cartridges	1	062948	48 Pack
Dionex OnGuard II Na Cartridges	2.5	062962	48 Pack

Dionex OnGuard II P Cartridges

- The Dionex OnGuard II P cartridge is recommended for removing the phenolic fraction of humic acids, tannic acids, lignins, anthocyanins, and azo dyes from samples prior to analysis by anion or cation exchange
- Cartridges contain polyvinylpyrrolidone (PVP) polymer with very high selectivity for phenolics, azo-containing compounds, aromatic carboxylic acids, aromatic aldehydes, and iodine as the triiodide complex

Key Specifications

- Functionality: polyvinylpyrrolidone
- Capacity ($\mu\text{eq}/\text{cartridge}$): 6.0 (1mL); 2.5mL format not available
- Solvents: 0-100% HPLC
- pH Stability: 1-10
- Mode: removal of phenols, azo dyes, humic acids by complexation

Dionex OnGuard II P Cartridges

Description	Volume (mL)	Cat. No.	Quantity
Dionex OnGuard II P Cartridges	1	057087	48 Pack

Dionex OnGuard II RP Cartridges

- The Dionex OnGuard II RP cartridge is recommended for removing hydrophobic substances such as aromatic dyes, some aromatic carboxylic acids, hydrocarbons, and surfactants from sample matrices
- Cartridges contain a macroporous divinylbenzene resin that has a very high selectivity for hydrophobic substances, especially unsaturated or aromatic organic substances

Key Specifications

- Functionality: polydivinylbenzene
- Capacity ($\mu\text{eq}/\text{cartridge}$): 0.3g resin (1mL); 0.75g resin (2.5mL)
- Solvents: 0-100% HPLC
- pH Stability: 0-14
- Mode: removal of surfactants, high-molecular weight carboxylic acids, aromatic dyes by adsorption

Dionex OnGuard II RP Cartridges

Description	Volume (mL)	Cat. No.	Quantity
Dionex OnGuard II RP Cartridges	1	057083	48 Pack
Dionex OnGuard II RP Cartridges	1	082760	12 Pack
Dionex OnGuard II RP Cartridges	2.5	057084	48 Pack

Dionex OnGuard Accessories

The Dionex OnGuard Sample Prep Station

- Enables simultaneous pretreatment of multiple samples with Dionex OnGuard sample pretreatment cartridges
- When used with 0.5mL PolyVials and a vacuum source, the Dionex OnGuard Sample Prep Station supports semi-automatic pretreatment of up to 12 samples
- Samples may also be treated manually using standard luer tip syringes
- Individual stopcock valves on each sample tube to allow control of individual flow rates
- The station will hold six 10mL volumetric flasks

Dionex OnGuard Accessories

Description	Cat. No.	Quantity
Dionex OnGuard Sample Prep Workstation	039599	1 Each
Dionex OnGuard Needle, 18 Gauge, 1.25/Luer	039996	1 Each
Dionex OnGuard Valve, Stopcock Luer	040896	1 Each

Dionex InGuard Cartridges

The Thermo Scientific™ Dionex™ InGuard™ line of sample pretreatment cartridges is designed to remove matrix interferences such as anions, cations, transition metals, or hydrophobic substances encountered in many IC and HPLC applications. The Dionex InGuard cartridge is installed inline between the autosampler and the IC injection valve facilitating immediate, automated sample pretreatment.

- Eliminates manual sample pretreatment steps
- Facilitates better separations
- Increases lifetimes of analytical columns
- Solves major matrix problems
- Achieves reproducible ppm-level determinations in concentrated matrices
- Convenient and easy to use



The Dionex InGuard cartridge is optimized for the best performance in matrix removal applications and can be used singly or in series. Designed to eliminate leaks and channeling, the cartridges use standard 10-32 fittings for easy installation into an IC or HPLC system. The unique sample distribution frit maximizes complete resin bed usage.

Note: Depending on the chemistry and samples treated, some cartridges can be regenerated.

Related Literature

For detailed product specifications and applications, see the following, available under Literature on www.thermoscientific.com/dionex

Product Datasheets

Dionex InGuard In-Line Sample Pretreatment Cartridges Data Sheet

Dionex InGuard Ag Cartridges

- Removes chloride, bromide, and iodide from concentrated sample matrices such as brines
- Styrene-based sulfonic acid resin in the silver form, the same material used in Dionex OnGuard II Ag cartridges
- For the removal of any residual silver ions, a Dionex InGuard H or Dionex InGuard Na cartridge should be placed after the Dionex InGuard Ag cartridge

Key Specifications

- Functionality: Cation-exchange, silver form
- Capacity: 5-5.5meq
- Solvents: 0-100% HPLC
- pH: 0-14
- Mode: Removal of halides by precipitation

Dionex InGuard Ag Cartridges

Description	Cat. No.	Quantity
Dionex InGuard Ag Cartridges	074038	4 Pack

Dionex InGuard H Cartridges

- The Dionex InGuard H cartridge is ideal for the removal of high levels of alkaline earth metals and transition metals from sample matrices. It is also used for the neutralization of highly alkaline samples such as sodium hydroxide or sodium carbonate
- Carbonate can be reduced to very low levels following this pH reduction by passing the sample through a Thermo Scientific™ Dionex™ CRD 200 Carbonate Removal Device
- Contains styrene-based, sulfonic acid resin in the hydronium form, the same as that used in Dionex OnGuard II H cartridges
- Resin is designed to have very high selectivity for polyvalent cations, such as calcium and transition metals

Key Specifications

- Functionality: Cation-exchange, hydronium form
- Capacity: 5-5.5meq
- Solvents: 0-100% HPLC
- pH: 0-14
- Mode: Removal of alkaline earth and transition metals; pH adjustment of basic samples

Dionex InGuard H Cartridges

Description	Cat. No.	Quantity
Dionex InGuard H Cartridges	074037	4 Pack

Dionex InGuard Na Cartridges

- The Dionex InGuard Na cartridge is used for the removal of high levels of alkaline earths and transition metals from sample matrices without acidifying the sample, ensuring good recovery of acid labile analytes such as nitrite
- Contains styrene-based, sulfonic acid in the sodium form, designed to have high selectivity for multivalent cations

Key Specifications

- Functionality: Cation-exchange, sodium form
- Capacity: 5-5.5meq
- Solvents: 0-100% HPLC
- pH: 0-14
- Mode: Removal of alkaline earth and transition metals

Dionex InGuard Na Cartridges

Description	Cat. No.	Quantity
Dionex InGuard Na Cartridges	074036	4 Pack

Dionex InGuard HRP Cartridges

- The Dionex InGuard HRP cartridge can be used to remove organic matrix material over a wide range of hydrophobicity, including fats from whole milk
- Contains a hydrophilic reversed-phase resin based on divinylbenzene. The material is water-wettable, thus 100% aqueous samples can be pretreated without disruption of the column bed

Key Specifications

- Functionality: Hydrophilic divinylbenzene
- Capacity: 2g
- Solvents: 0-100% HPLC
- pH: 0-14
- Mode: Adsorption, π - π bonding. Removal of hydrophobic species, azo-, and cyano-containing species

Dionex InGuard HRP Cartridges

Description	Cat. No.	Quantity
Dionex InGuard HRP Cartridges	074034	4 Pack

Dionex InGuard Na/HRP Cartridges

- The Dionex InGuard Na/HRP cartridge is designed to provide general purpose clean-up of samples, such as foods, for anion analysis
- Contains a blend of sulfonated resin in the sodium form and HRP resin to provide the dual functionality of removing both organic contaminants and cations, including metals, from a sample

Key Specifications

- Functionality: Dual Functionality
- Capacity: 50% Na/50% HRP
- Solvents: 0-100% HPLC
- pH: 0-14
- Mode: Ion-exchange (Na) and adsorption remove Ca^{2+} (Na^+) and lipids (HRP) from dairy products

Dionex InGuard Na/HRP Cartridges

Description	Cat. No.	Quantity
Dionex InGuard Na/HRP Cartridges	074035	4 Pack

Thermo Scientific Automated SPE

Dionex AutoTrace Systems

Provides fast and reliable automated solid phase extraction of organic pollutants from samples

The Thermo Scientific™ Dionex™ AutoTrace™ 280 instrument provides reliable automated SPE for analytical chemists determining organic pollutants in large volume liquid samples (10mL to 20L). Compared to liquid-liquid extraction, the Dionex AutoTrace 280 saves time, solvent, and labor, ensuring high reproducibility and productivity for analytical labs. The system uses powerful pumps (no check valves) and positive-pressure constant-flow technology to process the most difficult samples and can process up to six samples. Features include:

- SPE technology for liquid-liquid extraction: Reduces solvent usage and eliminates glassware.
- No technician involvement is required to maintain a liquid reservoir or to control the flow.
- Provides constant flow of liquids through SPE cartridges resulting in superior reproducibility.
- Closed systems with fan to vent solvent vapors: No hood required, conserves valuable hood space.
- The instrument can store 24 methods on board and the software used for editing methods is easy to use.
- USB cable connection to PC (methods to run can be chosen from those stored in the instrument memory).
- Multi-port switching valve ensures the systems reliability.



The Dionex AutoTrace 280 instrument provides automated solid-phase extraction with a choice of cartridge or disk formatting.

Dionex AutoTrace 280

Description	Cat. No.	Quantity
Dionex AutoTrace 280 Automated Large Volume SPE for 47mm Disks	071386	1 Each
Dionex AutoTrace 280 Automated Large Volume SPE for Cartridges with 1mL Plungers	072604	1 Each
Dionex AutoTrace 280 Automated Large Volume SPE for Cartridges with 3mL Plungers	072605	1 Each
Dionex AutoTrace 280 Automated Large Volume SPE for Cartridges with 6mL Plungers	071385	1 Each
Dionex AutoTrace 280 Automated Large Volume SPE for Cartridges with 6mL Glass Plungers	072606	1 Each



Dionex AutoTrace Accessories

The Dionex AutoTrace 280 parts and accessories are designed to be used specifically with the Dionex AutoTrace 280 solid-phase extraction system. Each part and accessory meets our strict quality standards.

SolEx SPE Cartridges

SolEx SPE Cartridges

Description	Bed Weight (mg)	Column Volume (mL)	Cat. No.	Quantity
SolEx C18	100	1	074623	100 Pack
SolEx C18	500	3	074412	30 Pack
SolEx C18	500	6	074417	30 Pack
SolEx C18	1,000	6	074410	30 Pack
SolEx C18 Unendcapped	1,000	6	074416	30 Pack
SolEx C18 Endcapped	1,000	6	075895	30 Pack
SolEx C18 Clean	1,500	6	075896	30 Pack
SolEx C8	100	1	074415	100 Pack
SolEx C8	500	3	074413	30 Pack
SolEx C8	1,000	6	074411	30 Pack
SolEx C8 Endcapped, Clean	500	6	075897	50 Pack
SolEx Silica	500	6	074589	50 Pack
SolEx Carbon-535	500	6	075898	30 Pack
SolEx Carbon-521	2,000	6	074590	20 Pack

Plunger Assemblies

Plunger Assemblies

Description	Cat. No.	Quantity
Plunger Assembly for 1mL Columns	071078	1 Each
Plunger Assembly for 3mL Columns	071079	1 Each
Plunger Assembly for 6mL Columns	071080	1 Each
Plunger Assembly for 6mL Glass Cartridges	071081	1 Each

Elution Racks

Elution Racks

Description	Cat. No.	Quantity
Elution Rack for 11mm GC vials	071068	1 Each
Elution Rack for 15mL Conical Tubes	071069	1 Each
Elution Rack for 16x100mL Test Tubes	071070	1 Each
Elution Rack for 17x60mm Vials	071071	1 Each
Elution Rack for 4mL Vials	071072	1 Each
15mL Conical Tubes	071056	Case of 12

Sample Rack

The Dionex AutoTrace 280 sample rack will hold six 60mL vials, 250mL bottles or 1L bottles. The rack is angled to ensure all the sample is retrieved from the sample vessels.

Sample Rack

Description	Cat. No.	Quantity
Sample Rack for 60, 250, and 1000mL Bottles for Dionex AutoTrace 280	071333	1 Each

Glass-Coated Solvent Bottles

Glass-coated solvent bottles are available in 1 and 2L sizes.

Glass-Coated Solvent Bottles

Description	Cat. No.	Quantity
Bottle, 1L, Coated Glass	045900	1 Each
Bottle, 2L, Coated Glass, GL45	045901	1 Each

Preventative Maintenance Kit

The preventative maintenance kit contains new tubing, seals and rotors.

Preventative Maintenance Kit

Description	Cat. No.	Quantity
Preventative Maintenance Kit	072598	1 Each

Thermo Scientific SPE Accessories

A comprehensive range of accessories to improve your sample preparation workflow

HyperSep Positive Pressure Manifold

HyperSep Positive Pressure Manifold

Description	Cat. No.	Quantity
Positive Pressure Manifold with 13mm Collection Rack	60104-236	1 Each
Positive Pressure Manifold with 16mm Collection Rack	60104-274	1 Each
Adaptor Plate for 1mL Cartridges for PP Manifold	60104-265	1 Each
Adaptor Plate for 3mL Cartridges for PP Manifold	60104-266	1 Each
Adaptor Plate for 6mL Cartridges for PP Manifold	60104-267	1 Each
Adaptor Plate for 10mL/15mL Cartridges for PP Manifold	60104-271	1 Each
Collection Rack for 13mm Tubes	60104-268	1 Each
Collection Rack for 16mm Tubes	60104-269	1 Each
Pre-Drilled Waste Container	60104-270	1 Each
Installation Kit	60104-272	1 Each
In Line Air Filter	60104-273	1 Each




HyperSep Glass Block Manifolds

16-Port Vacuum Manifold:

- Glass Block, Corian Manifold Lid, Cover Gasket, Vacuum Gauge and Valve Assembly, 16 Tips, Adjustable Collection Rack, Bulkhead Luer Fittings, 16 Plugs and Manifold Safety Tray

24-Port Vacuum Manifold:

- Glass Block, Corian Manifold Lid, Cover Gasket, Vacuum Gauge and Valve Assembly, 24 Tips, Adjustable Collection Rack, Bulkhead Luer Fittings, 24 Plugs and Manifold Safety Tray



Luer Lock Plugs



Retaining Clips for
Collection Rack



Bulkhead Luer Fittings



Collection Rack



Vacuum Gauge



Valve Assembly




Tips



Manifold Lid Legs

HyperSep Glass Block Manifolds

Description	Cat. No.	Quantity
16-port Vacuum Manifold	60104-232	1 Each
24-port Vacuum Manifold	60104-233	1 Each
Vacuum Pump, European Version	60104-241	1 Each
Vacuum Pump, North American Version	60104-243	1 Each
Replacement Parts		
Vacuum Gauge	60104-240	1 Each
Stopcocks for 16-port Vacuum Manifold	60104-242	16 Pack
Stopcocks for 24-port Vacuum Manifold	60104-244	24 Pack
Tips for Vacuum Manifold	60104-245	12 Pack
Lid for 16-port Glass Block Manifold	60104-262	1 Each
Lid for 24-port Glass Block Manifold	60104-248	1 Each
Gasket for 16-port Manifold	60104-249	1 Each
Gasket for 24-port Manifold	60104-250	1 Each
Collection Rack for 16-port Vacuum Manifold	60104-251	1 Each
Collection Rack for 24-port Vacuum Manifold	60104-252	1 Each
Glass Block for 16-port Vacuum Manifold	60104-253	1 Each
Glass Block for 24-port Vacuum Manifold	60104-254	1 Each
Retaining Clips for Collection Rack	60104-255	12 Pack
Bulkhead Luer Fittings	60104-256	12 Pack
Manifold Lid Legs	60104-257	4 Pack
Luer Lock Plugs	60104-258	12 Pack
Column Adaptors	60104-259	15 Pack
Manifold Safety Tray	60104-260	1 Each
Vacuum Gauge and Valve Assembly	60104-261	1 Each
Collection Rack Posts	60104-264	3 Pack



HyperSep Universal Vacuum Manifold

- Accommodates both SPE cartridges and 96 well plates
- System supplied with: Manifold, base/gauge, flask and stopper, tubing and spigots



Universal Manifold
Base/Gauge


HyperSep Universal Vacuum Manifold

Description	Cat. No.	Quantity
Universal Vacuum Manifold	60104-230	1 Each
Vacuum Pump, European Version	60104-241	1 Each
Vacuum Pump, North American Version	60104-243	1 Each
24 Well Extraction Plate	60104-284	1 Each
48 Well Extraction Plate	60104-289	1 Each
Plugs for 24-position Extraction Plate	60104-234	24 Pack
Plugs for 48-position Extraction Plate	60104-235	48 Pack

HyperSep-96 Well Plate Vacuum Manifold

Included with system:

- Base
- Lid
- Waste collection tray



HyperSep-96 well plate is not included and needs to be purchased separately.

HyperSep-96 Well Plate Manifolds

Description	Cat. No.	Quantity
HyperSep-96 Vacuum Manifold	60103-351	1 Each
Vacuum Pump, European Version	60104-241	1 Each
Vacuum Pump, North American Version	60104-243	1 Each

TurboFlow Columns

- For use with Thermo Scientific™ Transcend™ and Aria™ systems
- Perform turbulent flow chromatography
- Unique, patented column technology
- Automated online sample extraction
- Reduce costs, increase throughput and improve data quality



Thermo Scientific™ TurboFlow™ columns are used exclusively with the Transcend and Aria TLX Systems. Their patented design enables turbulent flow chromatography, which is valuable for online sample clean up by separating compounds of interest from larger molecular weight compounds or interferences in the sample matrix.

The compounds of interest are retained on the TurboFlow column while the unwanted matrix components flow quickly and completely to waste.

The retained compounds can then be transferred to an analytical column for further chromatographic separation prior to injection into a detector. All these steps are fully automated on the TLX systems.

Using TurboFlow technology, the Transcend and Aria TLX systems provide a powerful, high throughput sample preparation alternative to solid phase extraction (SPE), liquid-liquid extraction (LLE) and protein precipitation (PPT).

TurboFlow technology provides cleaner samples to the detector and thus, minimizes ion suppression, allows more reproducible data and reduces maintenance of the detector.

TurboFlow technology is ideal for automated online sample clean up of complex matrices; such as biological fluids and foods.

TurboFlow columns are available in a wide variety of chemistries to enable extraction of many different compounds of interest.

TurboFlow XL columns have an improved silica packing material and TurboFlow Cyclone columns are polymer based. Both columns are robust and have long lasting performance. All columns are available in both 1.0 and 0.5mm ID.

Transcend and Aria System TurboFlow Columns

Product Description	1x50mm	0.5x50mm
TurboFlow Cyclone Column	CH-952434	CH-953288
TurboFlow Cyclone-P Column	CH-952605	CH-953289
TurboFlow Cyclone MAX Column	CH-952979	CH-953286
TurboFlow Cyclone MCX Column	CH-952813	CH-953287
TurboFlow Cyclone MCX-2 Column	CH-953456	CH-953457
TurboFlow XL Column C18	CH-953244	CH-953280
TurboFlow XL Column C18-P	CH-953275	CH-953281
TurboFlow XL Column C2	CH-953279	CH-953285
TurboFlow XL Column C8	CH-953276	CH-953282
TurboFlow XL Column Fluoro	CH-953277	CH-953283
TurboFlow Column Phenyl*	CH-952178	CH-952820
TurboFlow Column PolarPlus™*	CH-952242	CH-952708

*Suggested operating pressure of 250 bar or below

For more information please visit www.thermoscientific.com/transcend

HyperSep Online SPE Products

Retain specific analytes in a sample matrix when used with an appropriate HPLC column

- Effective removal of contaminants such as proteins from samples
- Compatible with conventional HPLC systems
- Fast and effective clean-up and concentration of target compounds
- HyperSep Retain PEP: for retention of polar and non-polar analytes
- HyperSep Retain CX: for retention of basic and non-polar analytes
- HyperSep Retain AX: for retention of acidic and non-polar analytes
- HyperSep Hypercarb: for retention of extremely polar analytes



HyperSep Javelin Direct-Connect Online SPE Columns

ID (mm)	Length (mm)	Retain PEP	Retain CX	Retain AX	Hypercarb	Quantity
2.1	10	60310-201	60310-301	60310-401	60310-501	4 Pack
3.0	10	60310-202	60310-302	60310-402	60310-502	4 Pack


HyperSep UNIGUARD Direct-Connect Online SPE Cartridges

ID (mm)	Length (mm)	Retain PEP	Retain CX	Retain AX	Hypercarb	Quantity
2.1	10	60311-201	60311-301	60311-401	60311-501	4 Pack
3.0	10	60311-202	60311-302	60311-402	60311-502	4 Pack

HyperSep HPLC Columns for Online SPE

ID (mm)	Length (mm)	Retain PEP	Retain CX	Retain AX	Hypercarb	Quantity
2.1	20	60312-201	60312-301	60312-401	60312-501	1 Each
3.0	20	60312-202	60312-302	60312-402	60312-502	1 Each

This diagram shows the typical load and elution positions for the HyperSep Online SPE Setup



Dionex SolEx HRP Online Sample SPE Concentration Columns

Thermo Scientific™ Dionex™ SolEx™ Online Solid-phase extraction (SPE) columns are designed for fast, easy preconcentration of contaminants from water samples prior to HPLC analysis. The column is placed in the valve to automatically concentrate the analytes of interest. The online Dionex SolEx HRP columns are packed with a hydrophilic reversed-phase packing.

- Greater confidence in results due to less sample handling
- More reproducible results because of automated online use
- A wide range of compounds: from polar to hydrophobic analytes
- High capacity with good recoveries
- Wide pH compatibility (0-14)

Dionex SolEx HRP columns (2.1×20mm) are designed specifically for online HPLC use, and are packed with a divinylbenzene polymer with a hydrophilic bonded layer. This resin is used in the reversed-phase mode and demonstrates excellent retention properties on both polar and hydrophobic analytes. Even with 100% aqueous mobile phases, no de-wetting occurs. The Dionex SolEx HRP columns are available in both a cartridge format and an ultra high pressure compatible (UHPLC), rapid separation (RSLC) column format. The cartridge fits in a holder for easy cartridge exchange when necessary. The ultra high pressure, rapid separation column format can be used up to 800 bar.

Dionex SolEx Columns for Online HPLC Sample Preparation

Dionex SolEx Online Sample Preparation Columns	Cat No.
SolEx HRP cartridge (2.1x20mm), requires Holder V-3	074400
Cartridge Holder V-3	074403
SolEx RSLC HRP column (2.1x20mm)	075590

