

# HPLC Columns and Accessories

Solutions for liquid chromatography



**HAMILTON** 

# Innovation and excellence in chromatography

Hamilton Company has been developing and manufacturing pressure-stable polymeric high performance liquid chromatography (HPLC) columns for nearly 35 years. We are an established name in science whose products are found in most of the world's top chromatography labs. From columns to syringes to septa and more, Hamilton Company offers a full line of off-the-shelf and custom chromatography products for HPLC, gas chromatography (GC) and thin layer chromatography (TLC).

You can trust your results to Hamilton—The Measure of Excellence®



## Custom columns built to your exact need

See more on [page 48](#)



## HPLC syringes for every application

See more on [page 53](#)



## PRP-C18 columns: Excellent for extreme pH conditions

See more on [page 14](#)

# Table of Contents

- Hamilton: A Leader in HPLC .....3**
- Hamilton: Understanding Polymer Supports.....5**
- PRP-C18 Performance..... 8**
- HPLC Columns by USP Listing .....9**
- At-a-Glance Product Selection Charts..... 11**
- Reversed-Phase HPLC Columns..... 13**
  - PRP-C18 Columns ..... 14
  - PRP-1 Columns..... 17
  - PRP-3 Columns ..... 19
  - PRP-h5 Columns .....21
- Anion Exchange HPLC Columns .....23**
  - PRP-X100 and PRP-X110 Columns.....24
  - PRP-X500 Columns .....27
  - PRP-X600 Columns .....29
  - RCX-10 and RCX-30 Columns .....31
- Cation Exchange HPLC Columns.....33**
  - PRP-X200 Columns .....35
  - PRP-X400 Columns .....37
  - PRP-X800 Columns .....39
  - HC-40 Ca<sup>2+</sup> and HC-75 (H<sup>+</sup>, Ca<sup>2+</sup>, Pb<sup>2+</sup>) Columns..... 41
- Ion Exclusion HPLC Columns .....45**
  - PRP-X300 Columns .....45
- Custom HPLC Columns and Packing Materials...48**
- HPLC Guard Columns.....49**
- Bulk Polymer Resin.....50**
- Proper Column Care and Storage.....51**
- Chromatography Syringes .....53**
- About Hamilton Company .....54**

For more information on the full portfolio of Hamilton HPLC columns and other chromatography products, please visit [www.hamiltoncompany.com/HPLC](http://www.hamiltoncompany.com/HPLC) or refer to the back of this catalog for additional contact details.



# Hamilton: A Leader in HPLC

Founded on the invention of the Microliter™ syringe, Hamilton Company has been designing and innovating the industry's best syringes since 1953—not commercial, mass produced medical syringes, but syringes designed for precision sample measurement in chromatography instrumentation. Hamilton syringes have been trusted by chromatographers for nearly 60 years and for good reason—they are The Measure of Excellence.

The commitment to chromatography expanded into HPLC columns about 34 years ago. Hamilton was one of the first companies to understand the unique qualities of polymer-based columns and how the technology could advance the field of HPLC. Following significant research and advancements by its team of engineers and scientists, Hamilton's first line of pressure-stable polymeric HPLC columns was created, making it a pioneer in the development and manufacturing of polystyrene-divinylbenzene (PS-DVB) polymers for HPLC applications.

## Types of HPLC columns

Hamilton offers 17 different polymer-based HPLC columns for reversed-phase, anion exchange, cation exchange and ion exclusion separations, and two silica-based C8 and C18 columns for traditional reversed-phase separations.

The supports in a polymer-based column combine the inertness and pH stability of polymeric resins with the pressure stability and durability of silica-based materials. With Hamilton Polymeric Reversed-Phase (PRP™) HPLC columns and resins, the sample dictates the necessary separation conditions, not the limitations of the column.

Superior polymeric HPLC columns, resins and applications are a Hamilton specialty.

Nineteen different column packing materials are available for almost any challenging separation, including:

- ▶ Reversed-phase
- ▶ Anion exchange
- ▶ Cation exchange
- ▶ Ion exclusion

Specialty resins are available for a variety of difficult separations, including:

- ▶ Pharmaceuticals
- ▶ Herbicides
- ▶ Carbohydrates
- ▶ Proteins
- ▶ Alcohols



Need a custom method?

The Hamilton HPLC team is happy to help design your unique application. Give us a call to learn more.

## Bulk resins, column hardware and guard columns

Hamilton offers a variety of HPLC accessories to complement its column line.

- ▶ Hamilton HPLC column resins are available in bulk from 1 gram to kilogram quantities
- ▶ Column hardware in an array of specifications
  - ▶ Inside diameters: From 1.0–100 mm
  - ▶ Lengths: 30–305 mm
  - ▶ Hardware materials: Stainless steel and PEEK (polymer plastic)
- ▶ Guard columns are available to match the functionality of the analytical to preparative column sizes
  - ▶ Available dimensions: 3.0 x 8.0 mm, 2.0 x 20 mm, 4.6 x 20 mm
  - ▶ Available hardware materials: Stainless steel and PEEK

From autosampler syringes to manual injection and more, Hamilton has the chromatography syringe you need. View the full portfolio of HPLC, GC and TLC syringes at [www.hamiltoncompany.com/syringes](http://www.hamiltoncompany.com/syringes).



# Hamilton: Understanding Polymer Supports

What are polymer-based columns, and how are they different than silica?

Silica has its limitations. It requires functionalization for reversed-phase separations and is prone to chemical and pH degradation. Polymer resin manufacturing employs a unique process to produce, but it is this special science that makes them a truly viable and attractive alternative to silica-based materials. Polymer-based columns work comparably to silica-based columns, and, in many applications, perform even better. They can last for years, making them an easy and long-term cost savings alternative.

## Basic polymer structure

Polymer columns come in steel or PEEK hardware, just like silica-based columns, but the particles inside are a rigid polymer matrix rather than silica. Styrene (vinyl benzene) readily forms a polymer because the vinyl groups link together to form a chain. Cross-linking within the styrene groups occurs with the addition of divinyl benzene (DVB), which has a second vinyl group (i.e., meta or para to the first one). Cross-linking forms a much stronger and more rigid polymer. For HPLC applications, sufficient DVB is used to give a high density of cross-links thereby creating a more robust polymer support available for HPLC. This reaction is precisely controlled, allowing the formation of small spherical particles with a very narrow particle size distribution (within  $\pm 1 \mu\text{m}$ ).

## Advantages of polymers

Hamilton HPLC columns combine the best characteristics of silica-based and polymeric columns to arrive at a column that is highly inert and long-lasting. These characteristics are especially true for difficult analyses that require high pH (8–13), labile or reactive samples (e.g., irreversible adsorption), high aqueous purifications (80–100% water) and separations with ion pairing reagents.

## Alkaline pH stability

Mobile phase pH is a powerful tool in methods development, particularly for separation of neutral forms of amines or other organic bases under alkaline conditions. Although some recent C18 columns boast stability in alkaline pH, all silica-based supports experience measurable degradation at  $\text{pH} > 6$ , where column life is still considerably shorter than if used under more favorable conditions.

Polymeric columns, on the other hand, have genuine pH and chemical stability. The stationary phase stands up to prolonged exposure to concentrations as high as 1 M NaOH and  $\text{H}_2\text{SO}_4$ , with no measurable decrease in performance. Because the support does not strip, bleed, or dissolve at any pH, it can be expected to perform reliably and reproducibly throughout the extended life of the column, regardless of mobile phase conditions.

## Durable, long life polymers

Some polymers are prone to swelling in high organic solvents, rendering higher back pressures, but this is not the case with Hamilton's materials since a high degree of cross-linking prevents this from happening. Hamilton's PS-DVB supports are cross-linked to prevent shrinking or swelling, making them pressure stable up to 5,000 psi. Since the support is entirely polymeric with no silica to deteriorate, typical polymer column lifetime is approximately one year as compared to 3–4 months for an equivalent silica-based C8 or C18 column under routine or extreme method conditions.

## Long column life

The highly inert polymeric support resists chemical attack from organic solvents and aqueous buffers (0–100% aqueous or organic;  $\text{pH} 1\text{--}13$ ), effectively lengthening column life. If column performance deteriorates (e.g., peak broadening and a loss of symmetry), a regeneration protocol will usually return the chromatography to its original state. Even caustics can be used in extreme column fouling situations to reverse the adverse effects. This durability is especially important since it gives the chromatographer the ability to clean and regenerate a more expensive preparative column without the need to replace it.

## Wide application utility

The pH stability of Hamilton's polymer HPLC columns allows samples to be analyzed at basic pH (8–13). Some samples show a dramatic increase in their absorbance characteristics at basic pH. Altering the pH of the separation may not only increase a sample's detectability, it can also radically alter a compound's retention and the elution order of a sample. Changes in elution order can be used to determine sample purity and identity.



**Cleaning and more**

Polymer columns can be cleaned with 1.0 M sodium hydroxide or 100% organic solvent to remove strongly retained material and can be operated at much higher temperatures as compared to silica-based materials. They provide excellent sample recoveries due to the lack of acidic silanol groups associated with silica-based materials. Swelling in organic solvents such as THF or chloroform is negligible because of cross-linking.

PS-DVB resins are similar in retention characteristics to silica C18 but do have a slightly different selectivity in some cases.



# PRP-C18 Performance: Before and after harsh conditions


To demonstrate the rugged, pH-stability of Hamilton polymeric HPLC columns, a 10-component test mixture was run on a PRP-C18 column to obtain an initial separation. The column was then subjected to 200 column volumes of 1 M sodium hydroxide and then tested once again with the initial conditions to demonstrate that there was no measurable deterioration. After the sodium hydroxide flush, the column was then subjected to 200 column volumes of 1 M sulfuric acid and then tested once again with the initial conditions to demonstrate that there was still no measurable deterioration.

These tests demonstrate that polymer-based columns, unlike silica-based columns, can be freed of otherwise irreversibly bound contaminants under conditions that normally shorten column life.



**Hamilton is  
your partner.**


**From determining  
the correct column  
for your application  
to post-purchase  
support and  
troubleshooting,  
the HPLC team is  
standing by and  
ready to help.**



# HPLC Columns by USP Listing

This guide identifies all Hamilton HPLC columns according to their United States Pharmacopeial Convention (USP) listing.


- L17** Strong cation exchange resin consisting of sulfonated cross-linked styrene-divinylbenzene copolymer in the hydrogen form, 7 to 11  $\mu\text{m}$  in diameter
- ▶ PRP-X200..... Page 35
  - ▶ PRP-X300..... Page 45
  - ▶ HC-75 H<sup>+</sup>..... Page 41
- L19** Strong cation exchange resin consisting of sulfonated cross-linked styrene-divinylbenzene copolymer in the calcium form, about 9  $\mu\text{m}$  in diameter
- ▶ HC-75 Ca<sup>2+</sup>..... Page 41
  - ▶ HC-40 Ca<sup>2+</sup>..... Page 41
- L21** A rigid, spherical styrene-divinylbenzene copolymer, 3 to 10  $\mu\text{m}$  in diameter
- ▶ PRP-C18..... Page 14
  - ▶ PRP-1..... Page 17
  - ▶ PRP-3..... Page 19
  - ▶ PRP-h5..... Page 21
- L22** A cation exchange resin made of porous polystyrene gel with sulfonic acid groups, about 10  $\mu\text{m}$  in size
- ▶ PRP-X200..... Page 35
  - ▶ PRP-X300..... Page 45
- L23** An anion exchange resin made of porous polymethacrylate or polyacrylate gel with quaternary ammonium groups, about 10  $\mu\text{m}$  in size
- ▶ PRP-X500..... Page 27
- L34** Strong cation exchange resin consisting of sulfonated cross-linked styrene-divinylbenzene copolymer in the lead form, about 9  $\mu\text{m}$  in diameter
- ▶ HC-75 Pb<sup>2+</sup>..... Page 41
- L47** High capacity anion exchange microporous substrate, fully functionalized with trimethylamine groups, 8  $\mu\text{m}$  in diameter
- ▶ PRP-X100..... Page 24
  - ▶ PRP-X110..... Page 24
  - ▶ RCX™-10..... Page 31
  - ▶ RCX-30..... Page 31




# At-a-Glance Product Selection Charts

From molecular weight to separation type to product families, these simple, at-a-glance guides help determine the right Hamilton HPLC column for any need. For easy reference, each style of chromatography has a corresponding color that is also present in its product section within the catalog.

## Products by separation mechanism



## Products by molecular weight



# Reversed-Phase HPLC Columns

Hamilton reversed-phase HPLC columns combine the best characteristics of silica-based and polymeric columns to arrive at a product that is highly inert and long-lasting. Hamilton offers four polymeric packing materials for reversed-phase separations.

## Type Recommended Application(s)

- PRP-C18** Organic compounds: small molecules (< 2,000 mw), pharmaceuticals, steroids, organic halides, vitamins, amino acid analysis, herbicides
- PRP-1** Organic compounds: small molecules (< 2,000 mw), pharmaceuticals, steroids, nucleic acids, vitamins, herbicides
- PRP-3** Organic compounds: large molecules (> 2,000 mw), peptides, proteins, protein digests, protected and de-protected oligonucleotides, nucleic acids
- PRP-h5** Organic compounds: macromolecules (> 2,000 mw), pharmaceuticals, protein digests, tryptic digests

PS-DVB resins are similar in retention characteristics to silica C18s in that retention tends to increase with lipophilicity. However, subtle differences in the chemical interaction between the analyte and stationary phase can result in differential selectivity. In many cases, analytes that co-elute on a silica C18 can be resolved on a polymeric-based support. The PRP-1 consists of a 55% cross-linked PS-DVB bead containing 100 Å pores. The properties of this base material intrinsically lend itself to reversed-phase separations with no further surface modifications. The PRP-C18 uses the PRP-1 as the base support material with the addition of octadecyl to impart characteristics more closely related to a silica-based C18, giving slightly different selectivity than the PRP-1. To make the PRP-3, the PRP-1 is modified so that the base material contains 300 Å pores, which allows for the separation of larger molecules. The PRP-h5 utilizes the PRP-3 as its base with a pentafluorinated modification, making it more hydrophobic in nature.

In the early stages of reversed-phase chromatography, prototype columns were typically a silica bead functionalized with a C18 chain. These columns were revolutionary for their time but are being replaced by modern polymeric supports.

- ▶ PS-DVB supports are as retentive as silica C8 and C18 but offer alternate selectivity
- ▶ Stable over the full pH range (1–13)
- ▶ Compatible with virtually any aqueous, organic mobile phase
- ▶ Can be operated at temperatures well over 85°C
- ▶ Improved sample recovery compared to silica-based supports

## PRP-C18 Columns

### High efficiency separations at any pH

**Pore size: 100 Å**

**Material: C18-functionalized PS-DVB**

Mobile phase pH is a powerful tool in methods development, particularly for separation of neutral forms of amines or other organic bases under alkaline conditions. Although some recent C18 columns boast stability in alkaline pH, all silica-based supports experience measurable degradation at pH > 6, where column life is still considerably shorter than if used under more favorable conditions.

The PRP-C18, on the other hand, has genuine pH and chemical stability. The stationary phase stands up to prolonged exposure to concentrations as high as 1 M NaOH and H<sub>2</sub>SO<sub>4</sub>, with no measurable decrease in performance. Because the support does not strip, bleed, or dissolve at any pH, it therefore can be expected to perform reliably and reproducibly throughout the extended life of the column, regardless of mobile phase conditions.

### High pH applications

More than 70% of all pharmaceutical drug compounds are cationic solutes that carry a formal positive charge below pH 7. Separation of these and other organic bases has historically been problematic. Ionization has a dominating effect in reversed-phase chromatography that tends to dictate retention. Consequently, the elution window for a sample of ionized amines is narrow. The task is further complicated by secondary interactions that occur between positively charged solutes and residual silanols on the column stationary phase. These secondary mechanisms of retention are the principle source for anomalous chromatographic activity, such as poor peak shape, shifts in retention times and loss of efficiency that progressively worsen over the life of the column, as shown in the “Rapid Separation of Basic Drug Compounds on PRP-C18” chromatogram on page 15.

### Rapid elution

In modern drug discovery science, routine analytical chromatography should not be a bottleneck. As such, the trend is to increase productivity through the use of shorter columns packed with smaller particles and operated at elevated flow rates. The PRP-C18 is well suited for such use which is demonstrated in the chromatogram entitled “Separation of Common Organic Compounds on PRP-C18” on page 15.



### PRP-C18 stationary phase structure and applications

#### Applications:

Organic compounds: small molecules (< 2,000 mw), pharmaceuticals, steroids, halides, vitamins, amino acid analysis, herbicides


#### Examples of analytes that can be separated on PRP-C18 columns:

- ▶ Peptides
- ▶ DNA, RNA oligonucleotides, nucleotides
- ▶ Vitamins
- ▶ Steroids
- ▶ Herbicides
- ▶ Pharmaceutical compounds



### PRP-C18 application chromatograms


Separation of Common Organic Compounds on PRP-C18



**Column:** PRP-C18, 4.1 x 50 mm, 5 µm  
**Part number:** 79675  
**Mobile phase A:** Water + 0.2% H<sub>3</sub>PO<sub>4</sub>  
**Mobile phase B:** Acetonitrile + 0.2% H<sub>3</sub>PO<sub>4</sub>  
**Flow rate:** 2.5 mL/min  
**Gradient:** 2 to 99% B in 5 minutes  
**Temperature:** Ambient  
**Injection Volume:** 2 µL  
**Detection:** UV at 255 nm

- Compounds:**
1. Benzamide
  2. Nitromethane
  3. Nethyl 4-hydroxybenzoate
  4. N-ethyl 4-hydroxybenzoate
  5. N-propyl 4-hydroxybenzoate
  6. N-butyl 4-hydroxybenzoate
  7. Benzene
  8. Toluene
  9. Ethylbenzene
  10. Propylbenzene
  11. Pentylbenzene
  12. Hexylbenzene

Rapid Separation of Basic Drug Compounds on PRP-C18



**Column:** PRP-C18, 4.1 x 50 mm, 5 µm  
**Instrumentation:** Agilent 1100 quaternary pump with UV detector  
**Standards:**


1. Nicotine
2. Metoprolol
3. Quinine
4. Doxylamin
5. Dexmethorphan
6. Amitriptyline

**Mobile phase A:** 30 mM Diethylamine  
**Mobile phase B:** Mobile phase A + 95% ACN, 5% H<sub>2</sub>O  
**Gradient:** 10 to 100% B in 5 minutes  
**Flow rate:** 2 mL/min  
**Temperature:** Ambient  
**Injection volume:** 10 µL  
**Detection:** UV at 265 nm

### PRP-C18 Column Ordering Information

Dimensions	Particle Size		
	5 µm	10 µm	12-20 µm
2.1 x 50 mm	79672		
2.1 x 50 mm PEEK	79679		
2.1 x 150 mm	79673		
2.1 x 150 mm PEEK	79680		
2.1 x 250 mm	79674		
2.1 x 250 mm PEEK	79681		
4.6 x 50 mm	79675		
4.6 x 50 mm PEEK	79682		
4.6 x 150 mm	79676		
4.6 x 150 mm PEEK	79683		
4.6 x 250 mm	79677		
4.6 x 250 mm PEEK	79684		
21.2 x 250 mm			79678
Bulk Resin (1 gram)	79791	79792	79793

Learn more about PRP-C18 columns at [www.hamiltoncompany.com/PRPC18](http://www.hamiltoncompany.com/PRPC18).



View a keyword searchable index of applications possible with Hamilton HPLC columns at [www.hamiltoncompany.com/hplcapplicationindex](http://www.hamiltoncompany.com/hplcapplicationindex).



# PRP-1 Columns

## Superior sample recovery

**Pore Size:** 100 Å  
**Material:** PS-DVB

Sample recovery is vital to sample purification. Problems arise when labile samples become irreversibly bound to the silanol groups present on C8 and C18 HPLC columns. Since Hamilton polymers are made entirely of poly styrene-divinylbenzene, there are no silanol groups to cause sample loss. Recovery and quantitation of labile and reactive samples is enhanced. The purification of protected oligonucleotides demonstrates the enhanced recovery of polymer supports. While approximately 50–80% of an oligonucleotide is recovered on a C18 column, the equivalent PRP-1 column recovers 95% or greater of the same sample.


Unlike silica-based C8 or C18 columns, PRP-1 has no stationary phase coating. The integral reversed-phase characteristics of the PRP-1 column eliminate the need for special coating techniques. Since there is no stationary phase to hydrolyze, the column maintains its performance characteristics longer than many C8 or C18 columns.

### PRP-1 stationary phase structure and applications

**Applications:**  
Organic compounds: small molecule (< 2,000 mw), pharmaceuticals, steroids, nucleic acids, vitamins, herbicides

### Examples of analytes that can be separated on PRP-1 columns:

- ▶ Polycyclic aromatic hydrocarbons (PAH)
- ▶ Ionizable organic compounds
- ▶ Steroids
- ▶ Peptide fragments




Guard columns are an easy way to prolong an analytical column's life. Refer to page 49 for more information on how guard columns protect your investment.



### PRP-1 application chromatograms

Separation of Adenine and Guanine Nucleotides on PRP-1



**Column:** PRP-1, 5 µm, 2.1 x 150 mm  
**Part number:** 79366  
**Mobile phase A:** 100 mM Monopotassium phosphate (adjust pH to 7 with potassium hydroxide), 1 mM tetrabutylammonium phosphate, 2.5% methanol  
**Mobile phase B:** Mobile Phase A + 20% methanol  
**Flow rate:** 300 µL/min  
**Gradient:**  
Time (min) %B  
0–3 1  
10 15  
15 55  
16–19 95  
**Injection volume:** 5 µL  
**Sample concentration:** 0.02 mM  
**Temperature:** 50°C  
**Detection:** UV at 254 nm

**Compounds:**  
1. Guanosine monophosphate  
2. Guanosine diphosphate  
3. Adenosine monophosphate  
4. Guanosine triphosphate  
5. Adenosine diphosphate  
6. Adenosine triphosphate

### PRP-1 Column Ordering Information

Hardware Dimensions	Particle Size					
	5 µm	7 µm	10 µm	12–20 µm	30–50 µm	50–75 µm
1.0 x 50 mm		79755				
1.0 x 150 mm	79753					
2.1 x 100 mm	79790					
2.1 x 150 mm	79366					
4.1 x 50 mm	79443					
4.1 x 100 mm	79479					
4.1 x 150 mm	79444	79529	79425			
4.1 x 250 mm	79820	79422	79427			
4.6 x 100 mm PEEK	79558					
4.6 x 150 mm PEEK	79423		79351			
4.6 x 250 mm PEEK	79571	79380	79381			
7.0 x 100 mm			79495			
7.0 x 305 mm	79795		79426			
10.0 x 50 mm		79367				
10 x 100 mm	79355		79499			
10 x 250 mm		79531	79496			
21.2 x 75 mm	79154					
21.2 x 250 mm		79352	79478	79428		
30 x 250 mm				79229		
50 x 250 mm			79567	79493		
101.6 x 250 mm				79525		
101.6 x 250 mm Repack				79800		
Bulk Resin (1 Gram)	79578	79579	79580	79581	79582	79583



## PRP-3 Columns

Reversed-phase column optimized for separation of macromolecules such as DNA, RNA oligos, proteins, peptides and proteomics

**Pore size: 300 Å**  
**Material: PS-DVB**

Hamilton's PRP-3 is a polymeric reversed-phase HPLC column designed for the purification and isolation of proteins and peptides with very good recovery (> 90%). It is based off of the PRP-1 but utilizes a 300 Å pore size rather than the PRP-1's 100 Å pore size. The highly inert polymeric packing poly(styrene-divinylbenzene) enhances protein recovery because there are no silanol groups on the support to cause irreversible protein adsorption.

PRP-3 is a PS-DVB support that is pressure stable up to 5,000 psi and cross-linked to prevent shrinking or swelling when the mobile phase is changed. Chemically, proteins present solubility problems unlike many small molecules. Most proteins are hydrophobic on the inside, with highly charged exteriors. This often presents dissolution problems, particularly when pH is near the isoelectric point of the protein. The rugged chemical nature of the PRP-3 allows the protein chemist a much broader selection of agents for dissolution, including concentrated acids, aggressive chaotropes, as well as detergents.

In proteomics, the PRP-3 has excellent potential for single-column 2D HPLC. The orthogonal selectivity between low and high pH separations is often comparable to that achieved with two different column formats (SCX, RP), but with the added bonus of MS-compatible mobile phase.


### PRP-3 stationary phase structure and applications

#### Applications:

Organic compounds: large molecules (> 2,000 mw), peptides, proteins, protein digests, protected and de-protected oligonucleotides, nucleic acids


#### Examples of analytes that can be separated on PRP-3 columns:

- ▶ Globular proteins
- ▶ Albumins
- ▶ Antibody fragments
- ▶ Tryptic digests
- ▶ DNA
- ▶ RNA oligomers
- ▶ Synthetic high mw polymers



### PRP-3 application chromatograms

Five Proteins on PRP-3



Column: PRP-3, 4.1 x 150 mm, 10 µm

Part number: 79466

Standards:

1. Ribonuclease A
2. Cytochrome C
3. Lysozyme
4. Myoglobin
5. Ovalbumin

Mobile phase A: 0.01% TFA in water pH 2.0

Mobile phase B: 0.1% TFA in acetonitrile

Gradient: 25 to 50% B in 5 min. Hold 3 min.

Flow rate: 2 mL/min

Temperature: Ambient

Injection volume: 100 µL

Detection: UV at 215 nm

### PRP-3 Column Ordering Information

Hardware Dimensions	Particle Size	
	10 µm	12–20 µm
2.1 x 150 mm	79392	
4.1 x 150 mm	79466	
4.6 x 150 mm PEEK	79382	
4.6 x 250 mm PEEK	79574	
7.0 x 305 mm	79468	
10 x 250 mm	79526	
21.2 x 250 mm	79147	
21.2 x 100 mm		79186
21.2 x 250 mm		79469
Bulk Resin (1 gram)	79701	79702

A selection of column hardware sizes is available from analytical to semi-prep and preparative. Sample scale up is easy because the PRP-3 packing is consistent from analytical to preparative columns. This saves time and eliminates the need to redevelop separations on semi-prep or preparative columns. The short analytical (50 mm) column is well-suited for the gradient elution of high molecular weight proteins and the longer (150 mm) column is best for smaller proteins.

Custom HPLC columns are available! From dimensions to particle size to packing materials, Hamilton can build you exactly what you need. See page 48 for more information.



# PRP-h5 Columns

## Pentafluoro reversed-phase column for unique selectivity

**Pore size:** 300 Å

**Material:** Pentafluorinated PS-DVB

The PRP-h5 utilizes the PRP-3 as its base with a pentafluorinated modification, making it more hydrophobic in nature. The PRP-h5, with its functionality derived from a pentafluorinated polymer bead, delivers a selectivity difference from standard silica C18 stationary phases. This gives chromatographers the desired retention characteristics necessary for certain sample types. This difference is especially pronounced for halogenated solutes.

Unlike silica-based C8 or C18 columns, PRP-h5 has no stationary phase coating. Since there is no stationary phase to hydrolyze, the column maintains its performance characteristics longer than many silica-based C8 or C18 columns.


### PRP-h5 stationary phase structure and applications

#### Applications:


Organic compounds: macromolecules (> 2,000 mw), pharmaceuticals, protein digests, tryptic digests, proteomics

#### Examples of analytes that can be separated on PRP-h5 columns:

- ▶ Oligonucleotides
- ▶ Angiotensin
- ▶ apo-Transferrin
- ▶ Apomyoglobin (equine)
- ▶ Carbonic anhydrase
- ▶ Cytochrome C
- ▶ Myoglobin
- ▶ Ribonuclease A



### PRP-h5 application chromatograms



### PRP-h5 Column Ordering Information

Hardware Dimensions	Particle Size	
	5 µm	12-20 µm
2,1 x 100 mm	79270	
2.1 x 150 mm	79271	
4.6 x 50 mm	79261	
4.6 x 100 mm	79262	
4.6 x 150 mm	79272	
4.6 x 250 mm	79273	
10 x 100 mm	79263	
10 x 150 mm	79274	
Bulk Resin (1 gram)	79269	79280



Hamilton is your partner. From determining the correct column for your application to post-purchase support and troubleshooting, the HPLC team is standing by and ready to help.



# Anion Exchange HPLC Columns

Hamilton offers six polymeric packing materials for anion exchange separations.

Type	Recommended Application(s)
<b>PRP-X100</b> <b>PRP-X110</b>	Organic and inorganic anions, organic acids, organic and inorganic arsenic species
<b>PRP-X500</b>	Nucleic acids: single stranded/double stranded RNA and DNA Peptides and Proteins
<b>PRP-X600</b>	Adjustable exchange capacity by pH Nucleic acids: single stranded/double stranded RNA and DNA Peptides and Proteins
<b>RCX-10</b>	Carbohydrates, polysaccharides, sugar oligomers up to DP8
<b>RCX-30</b>	Mono and disaccharides

In anion exchange chromatography, the stationary bed has an ionically positive (+) charged surface while the sample ions are of negative (-) charge. This technique is used almost exclusively with ionic or ionizable samples. The stronger the negative charge on the sample, the stronger it will be attracted to the positive charge on the stationary phase, and thus the longer it will take to elute. Elution in ion chromatography is effected by mobile phase pH and ionic-strength, and, to a lesser extent, operation temperature. The ability to use the full pH range and elevated temperatures are distinct advantages compared to silica-based supports.

## PRP-X100 and PRP-X110 Columns

### Reliable separations of organic and inorganic anions

**Pore Size: 100 Å**

**Material: PS-DVB/Trimethyl ammonium exchanger**

**PRP-X100 Exchange Capacity: 0.19 meq/gm**

**PRP-X110 Exchange Capacity: 0.11 meq/gm**

Hamilton PRP-X100 and PRP-X110 are highly stable, inert materials. The PRP-X100 can be used with virtually any HPLC or ion chromatograph, including dedicated IC units. Technological advancements in modern polymer chemistry now deliver a more rugged column with exceptionally higher separation efficiencies than earlier predecessors. PRP-X100 and PRP-X110 columns are well suited for use in systems employing suppressed/non-suppressed conductivity, electrochemical, UV, and ICP-MS detection. Chromatographers currently using wet chemical or colorimetric methods will find ion chromatography greatly reduces sample pretreatment and improves the accuracy and precision of results.

PRP-X100 columns easily separate difficult anions such as cyanide, borate and silicate at high pH (11.5). The polymeric packing is stable from pH 1 to 13, so a single column can be used for the analysis of both common and difficult anions. The PRP-X100 is compatible with many different mobile phases for suppressed and non-suppressed conductivity and direct and indirect UV detection.

For high sensitivity, Hamilton PRP-X110 ion chromatography columns are used to separate ions at concentrations from less than 20 ppb to 20 ppm. The PRP-X110 has similar selectivity to the PRP-X100 but provides lower limits of detection as a result of its lower exchange capacity. The ion exchange capacity of a stationary phase plays a significant role in determining the concentrations of competing ions used in the mobile phase for elution. Lower capacity stationary phases generally require the use of weaker mobile phase eluent to affect elution due to the lower exchange capacity, thus improving the signal to noise. This is especially true when using conductivity detectors which do not function well with high salt eluents.

PRP-X110 columns can be used in the suppressed conductivity mode for determination of inorganic anions as required in EPA 300.0 Part A (fluoride, chloride, nitrite, bromide, nitrate, phosphate, sulfate). The PRP-X110 is a versatile column since it can be used with many different mobile phases such as carbonate, potassium hydroxide, benzoic acid, potassium hydrogen phthalate, etc. for suppressed and non-suppressed conductivity and direct and indirect UV detection.




## Stationary phase structure, applications and industries

### PRP-X100 and PRP-X110 applications:

Organic and inorganic anions, organic acids, organic and inorganic arsenic species

### PRP-X100 and PRP-X110 columns are used to monitor ions for a variety of industries including:

- ▶ **Pharmaceutical—**  
Lactate, acetate, chloride and phosphate in intravenous solution
- ▶ **Medical Research—**  
Monitoring analytes in bodily fluids of patients
- ▶ **Environmental—**  
Common anions in ground and river water, EPA 300.0 Part A
- ▶ **Food and Beverage—**  
Arsenic in food and water, phosphate in soft drinks and nitrates in food




### Examples of analytes that can be separated on PRP-X100 and PRP-X110 anion exchange columns:

- ▶ Halides (fluoride, chloride, etc.)
- ▶ Polarizable anions (perchlorate, thiocyanate)
- ▶ Organic acids, nucleotides, carboxylic acids (pyruvate, acetate, citrate, etc.)
- ▶ Organic and inorganic arsenic species

### PRP-X100 application chromatograms

Fast Separation of Common Anions on PRP-X100



Column: PRP-X100, 5  $\mu$ m, 2.1 x 50 mm

Part number: 79150

Mobile phase: 4 mM para-hydroxybenzoic acid, pH 8.5 / 4% acetonitrile

Flow rate: 0.8 mL/min

Gradient: Isocratic

Temperature: 80°C

Injection volume: 1  $\mu$ L


Sample concentration: 0.2 mg/mL

Detection: UV at 305 nm

#### Compounds:

1. Fluoride
2. Chloride
3. Nitrite
4. Bromide
5. Nitrate

Arsenic Speciation on PRP-X100



Column: Hamilton PRP-X100, 5  $\mu$ m, 4.6 x 250 mm

Part number: 79181

Flow rate: 1.0 mL/min

Mobile phase: 2 mM  $(\text{NH}_4)_2\text{CO}_3$  for 0–3 min  
40 mM  $(\text{NH}_4)_2\text{CO}_3$  for 3–14 min  
2 mM  $(\text{NH}_4)_2\text{CO}_3$  for 13–17 min

Injection volume: 50  $\mu$ L, 100  $\mu$ g/L of each standard

Detection: ICP-MS

#### Compounds:

1. Trivalent arsenic
2. Dimethyl arsenic
3. Monomethyl arsenic
4. Pentavalent arsenic
5. Monomethylthioarsonic acid

Borrowed (with permission) from P. Alava *et al.* Biomed. Chromatogr. 2012; 26: 524–533

### PRP-X100 and PRP-X110 Column Ordering Information

Hardware Dimensions	PRP-X100			PRP-X110	PRP-X110S
	Particle Size				
	5 $\mu$ m	10 $\mu$ m	12–20 $\mu$ m	7 $\mu$ m	7 $\mu$ m
2.1 x 100 mm PEEK				79743	
2.1 x 150 mm PEEK	79852				
2.1 x 150 mm		79421			
2.1 x 250 mm	79190	79346			
4.1 x 50 mm	79810	79365			
4.1 x 100 mm	79538	79439			
4.1 x 150 mm	79812	79434		79732	79733
4.1 x 250 mm		79433		79734	79735
4.6 x 150 mm PEEK	79174	79354		79738	
4.6 x 250 mm PEEK	79181	79455		79741	
100 Å 21.2 x 250 mm			79353		
Bulk Resin (1 Gram)	79584	79585	79586	79827	

### PRP-X100 Eluent Concentrate Ordering Information

Description	Part Number
Eluent Concentrate, PRP-X100 Anion Exchange (One 60 mL bottle)	79325
Eluent Concentrate, PRP-X100 Anion Exchange (Six 60 mL bottles)	79335

For a full list of all Hamilton HPLC products or for more detailed information, visit [www.hamiltoncompany.com/HPLC](http://www.hamiltoncompany.com/HPLC).

### The difference between PRP-X110 versus PRP-X110S columns

PRP-X110 columns are equilibrated with a 2 mM p-hydroxybenzoic acid pH 9.3 mobile phase and are ready for use with conductivity or indirect UV detection methods.

PRP-X110S columns are equilibrated with a 1.7 mM sodium bicarbonate, 1.8 mM sodium carbonate, 0.1 mM sodium thiocyanate mobile phase and are ready for use with suppressed conductivity detection methods.



## PRP-X500 Columns

### Fast separations and good sample recovery

**Pore Size:** Superficially porous

**Material:** Methacrylamido propyl trimethyl ammonium chloride (SAX)

**Exchange Capacity:** 1.6 meq/gm

PRP-X500 is a superficially porous polymeric anion exchange column designed for the separation, purification and isolation of proteins, peptides and DNA/RNA. The methacrylate polymeric coating of the PRP-X500 provides a more hydrophilic surface, preventing hydrophobic interaction sample losses typically seen on other commercially available protein HPLC columns.

The non-porous nature of the packing material improves mass-transfer, shortening run times and improving resolution. Both fast separations and good sample capacity are achievable with PRP-X500 columns.

A separation of four protein standards at 0.2 mg in less than three minutes is possible with a short analytical 50 x 4.6 mm HPLC column. Recovery of sample is excellent with PRP-X500 and the support's limited permeability prevents proteins from entering the pores and unfolding, which causes peak ghosting. The superficially porous properties shorten the diffusion path of the analyte, resulting in sharp sample bands.

### PRP-X500 stationary phase structure and applications


#### Applications:

Proteins/Peptides

Single Stranded/Double Stranded RNA/DNA


#### Examples of analytes that can be separated on PRP-X500 columns:

- ▶ Myoglobin
- ▶ Bovine serum albumin
- ▶ Conalbumin
- ▶ Ovalbumin



### PRP-X500 application chromatograms

Myoglobin, Conalbumin and Dog Albumin on PRP-X500



Column: PRP-X500, 4.6 x 50 mm, 5  $\mu$ m  
 Part number: 79474  
 Mobile phase A: 10 mM Tris pH 9.0  
 Mobile phase B: 70 mM Tris pH 9.0, 0.5 N Sodium Chloride  
 Flow rate: 2 mL/min  
 Gradient: 0 to 50% B in 2.5 minutes. Hold for 2.5 min.  
 Temperature: Ambient  
 Injection volume: 30  $\mu$ L  
 Detection: UV at 254 nm

Compounds:  
 1. Myoglobin 7  $\mu$ g  
 2. Conalbumin 7  $\mu$ g  
 3. Dog albumin 77  $\mu$ g

### PRP-X500 Column Ordering Information

Dimensions	Particle Size		
	7 $\mu$ m	12–20 $\mu$ m	30–50 $\mu$ m
4.6 x 150 mm, PEEK	79573		
Bulk Resin (1 gram)	79594	79595	79596

Guard columns are an easy way to prolong an analytical column's life. Refer to page 49 for more information on how guard columns protect your investment.



## PRP-X600 Columns

### Separation of nucleic acids

**Pore Size: Superficially porous**

**Material: Poly (dimethylamidopropylmethacrylamide)**

**Exchange Capacity: 1.6 meq/gm**

PRP-X600 is a superficially porous weak-base anion exchange support that separates DNA oligomers according to negative charge. The unique porosity provides fast separation with better sample capacity than non-porous supports. The superficially porous properties shorten the diffusion path of the analyte, resulting in sharp sample bands. The PRP-X600, a hydrophilic methacrylate-based polymer, improves sample recovery due to minimized hydrophobic interactions.

Because the PRP-X600 is a weak anion exchange (WAX) resin, the exchange capacity of the resin is pH dependent. Lowering the pH will reduce the binding of proteins, reducing the run time for a complete separation.

Change the mobile phase composition to alter the retention of DNA oligomers. Rapid gradient changes typically lower column efficiency; however, biomolecules run in this fashion on the PRP-X500 show very favorable separation efficiency with much shorter run times.


#### PRP-X600 stationary phase structure and applications

##### Applications:

Nucleic acids such as single stranded/double stranded RNA and DNA Peptides and proteins


##### Example of analytes that can be separated on the PRP-X600:

- ▶ Synthetic RNA, DNA oligonucleotides
- ▶ Proteins and peptides
- ▶ Ovalbumin
- ▶ DNA fragments
- ▶ Oligonucleotides



#### PRP-X600 application chromatograms

Oligodeoxycytidylate (dc) 12-18 on PRP-X600



**Column:** PRP-X600, 7  $\mu$ m, 4.6 x 50 mm  
**Part number:** 79360  
**Mobile phase A:** 85/15 100 mM TRIS, pH 8.0/acetonitrile  
**Mobile phase B:** 85/15 100 mM TRIS, pH 8.0, 2.5 M lithium chloride/acetonitrile  
**Flow rate:** 2.0 mL/min  
**Gradient:** 0 to 40% B in 40 minutes  
**Temperature:** Ambient  
**Injection volume:** 10  $\mu$ L  
**Sample concentration:** 300  $\mu$ g/mL  
**Detection:** UV at 260 nm

##### Compounds (from oligodeoxycytidylate (dC)<sub>12-18</sub>)

1. dC12
2. dC13
3. dC14
4. dC15
5. dC16
6. dC17
7. dC18

#### PRP-X600 Column Ordering Information

Dimensions	Particle Size		
	7 $\mu$ m	12–20 $\mu$ m	30–50 $\mu$ m
4.6 x 50 mm PEEK	79360		
4.6 x 250 mm PEEK	79189		
Bulk Resin (1 Gram)	79597	79598	79599



## RCX-10 and RCX-30 Columns

Designed for the isocratic or gradient separation of carbohydrates

**Pore Size:** 100 Å

**Material:** PS-DVB/Trimethyl ammonium exchanger

**RCX-10 Exchange Capacity:** 0.35 meq/gm

**RCX-10 Exchange Capacity:** 1.0 meq/gm

The Hamilton RCX-10 and RCX-30 carbohydrate analysis columns are designed for the isocratic or gradient separation of carbohydrates. The exchange capacity of the RCX-10/RCX-30 is greater than that of the PRP-X100, leading to characteristics better suited for the separation of carbohydrates. Simple samples with two or three carbohydrates can be quickly separated isocratically, while more complex samples require gradient elution to fully resolve all the analytes of interest. When an isocratic method is used with a conductivity, refractive index, ultraviolet, or pulsed amperometric detector (PAD), mono and disaccharides such as glucose, fructose, sucrose and lactose can be quickly determined.

To utilize the full potential of the RCX-10 column (e.g., gradient separations) a Pulsed Amperometric Detector (PAD) is recommended. The PAD allows utilization of either gradient or isocratic elution for the separation of carbohydrates in foods or food products.

A typical mobile phase is sodium hydroxide and sodium acetate. When the concentration of these mobile phases is varied, a variety of samples can be separated. Separation of carbohydrates with the RCX-10 or RCX-30 at basic pH is possible since each carbohydrate carries a different negative charge at basic pH.

RCX-30 carbohydrate analysis columns provide longer sample retention than RCX-10 columns and better resolution of complex samples like the six constituent monosaccharides of glycoconjugates. It is the increased exchange capacity of the RCX-30 that gives it these characteristics as compared to the RCX-10.


### RCX-10 and RCX-30 stationary phase structure and applications

#### Applications:



Carbohydrates, polysaccharides, sugar oligomers up to DP8 mono and disaccharides

#### Example of analytes that can be separated on RCX-10 and RCX-30 columns:

- ▶ Arabinose
- ▶ Galactose
- ▶ Lactose
- ▶ Maltose
- ▶ Sucrose



### RCX-10 and RCX-30 application chromatograms



### RCX-10 and RCX-30 Column Ordering Information

Dimensions	RCX-10		RCX-30	
	Particle Size			
	7 µm	12–20 µm	7 µm	12–20 µm
4.1 x 250 mm	79440		79803	
4.6 x 150 mm, PEEK			79370	
4.6 x 250 mm, PEEK	79388		79877	
Bulk Resin (1 Gram)	79703	79704	79705	79706

View a keyword searchable index of applications possible with Hamilton HPLC columns at [www.hamiltoncompany.com/hplcapplicationindex](http://www.hamiltoncompany.com/hplcapplicationindex).

# Cation Exchange HPLC Columns

Hamilton offers seven polymeric packing materials for cation exchange separations.

Type	Recommended Application(s)
<b>PRP-X200</b>	Inorganic and organic cations using conductivity or UV detection, alkali and alkaline earth metals. Separate mono or divalent cations depending on mobile phase conditions from 20 ppb to 200 ppm.
<b>PRP-X400</b>	Glyphosate and its metabolite in drinking water Organic and inorganic cations using conductivity or UV detection
<b>PRP-X800</b>	Transition metals or mono and divalent cations in the same run
<b>HC-40</b>	Oligo saccharides up to DP8
<b>HC-75 Ca<sup>2+</sup></b>	Mono and disaccharides
<b>HC-75 H<sup>+</sup></b>	Organic acids, sugars and alcohols
<b>HC-75 Pb<sup>2+</sup></b>	Sugar alcohols and plant cell wall hydrolysates

In cation exchange chromatography, the stationary bed has an ionically negative (-) charged surface while the sample ions are of positive (+) charge. This technique is used almost exclusively with ionic or ionizable samples. The stronger the positive (+) charge on the sample, the stronger it will be attracted to the negative charge on the stationary phase, and thus the longer it will take to elute. The mobile phase is an aqueous buffer, where both pH and ionic strength are used to control elution time. Ion chromatography can employ harsh conditions requiring mobile phases that are at very high pH limits (> 11). Temperatures well above the normal operating conditions where silica materials fail can also be used.



## PRP-X200 Columns

### High resolution separation of alkali and alkaline earth metals

**Pore Size:** 100 Å

**Material:** PS-DVB/Sulfonic acid exchanger

**Exchange Capacity:** 35 µeq/gm

Hamilton PRP-X200 cation exchange HPLC columns are designed for rapid, high resolution separation of alkali and alkaline earth metals. The alkali metals and ammonium are completely resolved in less than five minutes, and the alkaline earth cations separate in under four minutes. Since the mobile phase conditions are different and unique for each of the groups of cations, interferences between these groups are eliminated.

The resolution between any two ions in the alkali metal series can be increased or decreased by changing the concentration of methanol in the mobile phase. This unique feature of the PRP-X200 mobile phase/stationary phase interaction allows the chromatographer to focus on the particular ion of interest and reduce possible interference from other ions in the sample.

The stationary phase is a sulfonated poly (styrene-divinylbenzene), so it is stable to all concentrations of organic modifiers as well as strong acids and bases. The rigid spherical polymer phase allows operation up to 5,000 psi.

The PRP-X200 columns are designed for use in all of today's ion chromatographic equipment. Since exchange capacity is low, background signal is low. This allows low detection limits at the highest conductivity detector settings.


### PRP-X200 stationary phase structure and applications

#### Applications:

Inorganic and organic cations using conductivity or UV detection, alkali and alkaline earth metals. Separate mono or divalent cations depending on mobile phase conditions from 20 ppb to 200 ppm.


#### Examples of analytes that can be separated on PRP-X200 columns:

- ▶ Calcium
- ▶ Potassium
- ▶ Cesium
- ▶ Sodium



### PRP-X200 application chromatograms

*Monovalent Cations on PRP-X200*



**Column:** PRP-X200, 5 µm, 4.1 mm x 150 mm  
**Part number:** 79441  
**Mobile phase:** (2.3:1) 4 mM Nitric acid:Methanol  
**Flow rate:** 2 mL/min  
**Gradient:** Isocratic  
**Temperature:** Ambient  
**Injection volume:** 100 µL  
**Detection:** Conductivity

**Compounds:**  
 1. Lithium  
 2. Sodium  
 3. Ammonium  
 4. Potassium  
 5. Cesium



**Need a custom method? The Hamilton HPLC team is happy to help design your unique application. Give us a call to learn more.**

The 150 mm column is ideal for rapid analysis, while the 250 mm column is recommended for higher resolution or for analysis of minor components in the presence of major interferences.

#### PRP-X200 Column Ordering Information

Dimensions	Particle Size	
	10 µm	12–20 µm
2.1 x 150 mm	79394	
4.1 x 100 mm	79363	
4.1 x 150 mm	79441	
4.1 x 250 mm	79442	
4.6 x 150 mm PEEK	79384	
4.6 x 250 mm PEEK	79357	
Bulk Resin (1 Gram)	79587	79588

## PRP-X400 Columns

### Excellent separation of glyphosate and its metabolite in drinking water

**Pore Size:** 100 Å

**Material:** PS-DVB/Sulfonic acid exchanger

**Exchange Capacity:** 2.5 meq/gm

The PRP-X400 column provides a fast separation for glyphosate and its metabolites. The exchange capacity of the PRP-X400 is greater than that of the PRP-X200, leading to characteristics better suited for the separation of glyphosate. It also performs well in other separations, such as inositol and sugar alcohols. The column does not have to be heated to 65°C and operates well at room temperature, so a column heater is not necessary for this method. PRP-X400 columns do not require the use of methanol in the mobile phase, and they cost much less than other glyphosate columns.

The PRP-X400 is a 7 µm poly(styrene-divinylbenzene) sulfonated cation exchange support (2.5 meq/gm) column. It separates glyphosate and aminomethylphosphonic acid according to charge in less than 10 minutes. This separation requires post-column oxidation and derivatization.

Post column reaction (oxidation) with calcium hypochlorite followed by derivatization with o-phthalaldehyde solution provides sensitive (6 ppb or lower) and selective (primary and secondary) amine detection.


### PRP-X400 stationary phase structure and applications

#### Applications:

Glyphosate and its metabolite in drinking water. The PRP-X400 provides unique hydrophilic interaction separations. Cations, inorganic and organic using conductivity or UV detection.


#### Examples of analytes that can be separated on PRP-X400 columns:

- ▶ Glyphosate
- ▶ Xylitol
- ▶ Maltose
- ▶ Mannitol



### PRP-X400 application chromatograms

Glyphosate on PRP-X400



**Column:** PRP-X400, 5 µm, 4.1 mm x 250 mm  
**Part number:** 79473  
**Mobile phase:** 0.005 M Monobasic potassium phosphate  
**Flow rate:** 0.5 mL/min  
**Gradient:** Isocratic  
**Temperature:** Ambient  
**Injection volume:** 200 µL  
**Detection:** Excitation wavelength—338 nm, Emission wavelength—455 nm

**Compounds:**  
 1. Glyphosate  
 2. Aminomethylphosphonic acid

### PRP-X400 Column Ordering Information

Dimensions	Particle Size		
	7 µm	12–20 µm	30–50 µm
2.1 x 250 mm	79398		
4.1 x 150 mm	79717		
4.1 x 250 mm	79473		
4.6 x 250 mm PEEK	79387		
Bulk Resin (1 Gram)	79591	79592	79593

For full details on how to use the PRP-X400 for glyphosate analysis, please visit our website at [www.hamiltoncompany.com/HPLC](http://www.hamiltoncompany.com/HPLC).



## PRP-X800 Columns

Separates transition metals or mono and divalent cations in the same run

**Pore Size:** 100 Å

**Material:** PS-DVB/Itaconate exchanger (WCX)

**Exchange Capacity:** 1.6 meq/gm

The PRP-X800 is a polymeric cation exchange column functionalized with itaconic acid that performs the isocratic separation of mono and divalent cations such as lithium, sodium, ammonium, potassium, magnesium and calcium. The column offers excellent durability, is stable to any concentration organic solvent, and enables dynamic control of exchange capacity. Detection is via conductivity or indirect UV, depending on the mobile phase.


### PRP-X800 stationary phase structure and applications

#### Applications:

Mono and divalent transition metals in the same run. Transition metals (e.g., manganese, zinc, cobalt and cadmium) are also resolved on the column using an ethylenediamine/tartaric acid mobile phase and conductivity detection.


#### Examples of analytes that can be separated on PRP-X800 columns:

- ▶ Mono and divalent metal cations (e.g., sodium, potassium, calcium)
- ▶ Transitions metals (e.g., iron, manganese, nickel, copper, zinc)



### PRP-X800 application chromatograms

Mono and Divalent Cations on PRP-X800



**Column:** PRP-X800, 5 µm, 4.1 mm x 250 mm  
**Part number:** 79828  
**Mobile phase:** 2 mM Cupric Sulfate  
**Flow rate:** 0.8 mL/min  
**Gradient:** Isocratic  
**Temperature:** Ambient  
**Injection volume:** 10 µL  
**Sample concentration:** All compounds are 5 ppm  
**Detection:** Indirect UV at 220 nM

#### Compounds:

1. Lithium
2. Sodium
3. Ammonium
4. Potassium
5. Magnesium
6. Calcium

The PRP-X800 is available in PEEK hardware as a custom order. Please see page 48 for more information.

### PRP-X800 Column Ordering Information

Dimensions	Particle Size: 7 µm
4.1 x 150 mm	79855
4.1 x 250 mm	79828

From autosampler syringes to manual injection and more, Hamilton has the chromatography syringe you need. View the full portfolio of HPLC, GC and TLC syringes at [www.hamiltoncompany.com/syringes](http://www.hamiltoncompany.com/syringes).

## HC-40 $\text{Ca}^{2+}$ and HC-75 ( $\text{H}^+$ , $\text{Ca}^{2+}$ , $\text{Pb}^{2+}$ ) Columns

### Separate compounds through size exclusion and ligand exchange

**Pore Size:** 100 Å

**Material:** PS-DVB/ Sulfonic acid exchanger

**Exchange Capacity:** 5 meq/gm

The HC-40 and HC-75 group of columns consist of four different packing types, each with a different retention characteristic and application. They are:

- ▶ HC-40—Oligo saccharides up to DP8 (e.g., corn syrup, high conversion corn syrup, beer)
- ▶ HC-75 Lead Form—Sugar alcohols and plant cell wall hydrolysates
- ▶ HC-75 Calcium Form—Separates mono and disaccharides (e.g., corn syrup, chewing gum sweeteners, milk product sugars, glycols and polyols, high fructose corn syrup, juices, oligosaccharides)
- ▶ HC-75 Hydrogen Form—Organic acids, sugars and alcohols

HC-40 and HC-75 columns separate compounds through size exclusion and ligand exchange. The 4% cross-linked HC-40 uses size exclusion as the primary mechanism of separation, while ligand exchange dominates in the more highly cross-linked HC-75. The different forms of the HC-75 (Hydrogen, Calcium and Lead) each provide a unique selectivity for separating varying types of charged analytes based on electronegativity toward the counterion. The higher carbohydrate oligomers elute first while the smaller di- and monosaccharides elute later.

The HC-75 column provides a slightly faster 14-minute separation up to DP 5, and the HC-40 column provides a much better separation of the oligomers up to DP 8 in 16 minutes. Because carbohydrates do not contain a chromophore, UV detection cannot be used without derivatization. The recommended detection method is refractive index. The control of carbohydrate retention lies in the selection of the correct column.

The HC-40 and HC-75 columns use water as a mobile phase (gradients and salts are not required), which simplifies eluent preparation and minimizes cost. This mobile phase characteristic also lends these columns to detection techniques such as evaporative light scattering detection and mass spectrometry. A very versatile column family, the HC-75 group can be used with up to 40% acetonitrile and can be regenerated to help restore performance.


### HC-40 and HC-75 stationary phase structure and applications

#### Applications:

Carbohydrates, sugar oligomers up to DP8. Mono and disaccharides, organic acids, sugars, and sugar alcohols.

#### Examples of analytes that can be separated on HC-40 and HC-75 columns:

- ▶ Ethanol
- ▶ Maltotriose
- ▶ Glucose
- ▶ Fructose
- ▶ Arabinose
- ▶ Sorbitol
- ▶ Acetic acid




$\text{Pb}^{2+}$ ,  $2\text{H}^+$  or  $\text{Ca}^{2+}$



HC-75 application chromatograms


Chewing Gum Sugars on HC-75 Ca<sup>2+</sup>



**Column:** HC-75 Calcium Form, 5 μm, 7.8 mm x 305 mm  
**Part number:** 79436  
**Mobile phase:** Deionized water.  
**Flow rate:** 1.2 mL/min  
**Gradient:** Isocratic  
**Temperature:** 90°C  
**Injection volume:** 2 μL  
**Detection:** Refractive index

**Compounds:**  
 1. Sucrose  
 2. Glucose  
 3. Fructose  
 4. Mannitol  
 5. Sorbitol

Organic Acids by USP L17 on HC-75 H<sup>+</sup>




**Column:** HC-75 Hydrogen Form, 5 μm, 4.1 mm x 250 mm  
**Part number:** 79476  
**Mobile phase:** 0.01 N sulfuric  
**Flow rate:** 0.35 mL/min  
**Gradient:** Isocratic  
**Temperature:** 60°C  
**Injection volume:** 10 μL

**Sample concentration:** All compounds are 250 ppm  
**Detection:** UV at 210 nm

**Compounds:**  
 1. Citric acid  
 2. Lactic acid  
 3. Acetic acid

Sugar Standards on HC-75 Pb<sup>2+</sup>



**Column:** HC-75 Lead Form, 5 μm, 7.8 mm x 305 mm  
**Part number:** 79438  
**Mobile phase:** Deionized water  
**Flow rate:** 0.6 mL/min  
**Gradient:** Isocratic  
**Temperature:** 80°C  
**Injection volume:** 200 μL  
**Sample concentration:** All compounds are 2.5 mg/mL  
**Detection:** Refractive index

**Compounds:**  
 1. Sucrose  
 2. Glucose  
 3. Fructose  
 4. Inositol  
 5. Sorbitol

HC-75 and HC-40 Column Ordering Information

	HC-75	HC-40
<b>Dimensions</b>	<b>Particle Size</b>	
	9 μm	9 μm
(H <sup>+</sup> ) 4.1 x 250 mm	79476	
(H <sup>+</sup> ) 7.8 x 100 mm	79547	
(H <sup>+</sup> ) 7.8 x 305 mm	79544	
(H <sup>+</sup> ) Bulk Resin (1 Gram)	79711	
(Ca <sup>2+</sup> ) 4.1 x 250 mm	79431	
(Ca <sup>2+</sup> ) 7.8 x 305 mm	79436	
(Ca <sup>2+</sup> ) Bulk Resin (1 Gram)	79709	79707
(Pb <sup>2+</sup> ) 7.8 x 100 mm	79240	
(Pb <sup>2+</sup> ) 7.8 x 305 mm	79438	
(Pb <sup>2+</sup> ) Bulk Resin (1 Gram)	79712	



# Ion Exclusion HPLC Columns

Hamilton offers one polymeric packing material for ion exclusion separations.

Ion exclusion chromatography is an alternative to ion exchange chromatography in which ionized samples are excluded from the pores of the support and elute first, while the weakly ionized and nonionic compounds elute later. Mixtures of weak acids, like those in fruits and milk products, are frequently not very well separated by pure ion-exchange methods, nor in the reversed-phase mode.

## PRP-X300 Columns

### Ion exclusion for organic acids and alcohols

**Pore Size:** 100 Å

**Support Material:** PS-DVB/Sulfonic acid

Hamilton PRP-X300 columns offer an easy, rapid way to separate closely related alcohols and organic acids. The sulfonated poly(styrene-divinylbenzene) support separates samples via a mixed mode mechanism. Separation on the PRP-X300 is accomplished by three modes:

1. **Hydrogen Bonding**—The attraction and retention of sample compounds by the negatively charged sulfonate group.
2. **Reversed-Phase**—The interaction and retention of the sample compounds by the non-polar polymeric support.
3. **Ion Exclusion**—The process in which ionized samples are excluded from the pores of the support and elute first, while the weakly ionized and nonionic compounds elute later.

A wide variety of samples can be analyzed with PRP-X300 columns because their selectivity can be altered by changing the pH of the buffer or adding an organic modifier (e.g., methanol, acetonitrile). The support's stability to organic solvents makes it possible to analyze samples that are too highly retained on conventional ion exclusion supports.

Most separations on PRP-X300 columns are completed within five minutes. The high performance packing ensures narrow peaks under isocratic conditions. The use of isocratic conditions allows samples to be analyzed one after another without waiting for column re-equilibration. Complex samples can be run under gradient conditions if isocratic methods are insufficient. Most samples require only minimal preparation before injection, increasing sample throughput in the lab.


## PRP-X300 stationary phase structure and applications

### Applications:


Organic acids and alcohols

### Examples of analytes that can be separated on PRP-X300 columns:

- ▶ Acetic acid
- ▶ Acrylamide
- ▶ Citric acid
- ▶ Oxalacetic acid
- ▶ Ethanol
- ▶ Propanol




## PRP-X300 application chromatogram



## PRP-X300 Column Ordering Information

Dimensions	Particle Size	
	7 μm	12–20 μm
4.6 x 150 mm PEEK	79475	
4.1 x 250 mm	79465	
4.1 x 150 mm	79464	
Bulk Resin (1 Gram)	79589	79590

View an index of applications possible with Hamilton HPLC columns at [www.hamiltoncompany.com/hplcapplicationindex](http://www.hamiltoncompany.com/hplcapplicationindex).



# Custom HPLC Columns and Packing Materials

Hamilton offers an extensive line of off-the-shelf HPLC columns with different packing material chemistries, varying particle sizes and column hardware dimensions in two hardware types—stainless steel and polyether ether ketone (PEEK) plastic. However, even with all the standard options available, customers often require a custom built column where they may need one of our stationary phases packed into a column hardware combination that is not a standard product.

Custom HPLC columns offer customers the ability to have any of our stationary phases packed into most any hardware format to suit their application requirement.

## Ordering custom HPLC products

Customers can specify the following information either through the online shopping cart or over the phone with a customer service representative:

- ▶ Packing material and particle size (e.g., PRP-1, 5  $\mu\text{m}$ )
- ▶ Hardware dimensions and material (e.g., 4.6 x 150 mm, PEEK)

## Custom HPLC Column Ordering Information

Custom Column Category	Part Number
Analytical HPLC Column, Custom Order (1.0, 2.1, 4.1, 4.6 mm ID)	79641
Semi-preparative HPLC Column, Custom Order (7.0, 7.8, 10 mm ID)	79642
Preparative HPLC Column #1, Custom Order (21.2, 30 mm ID)	79643
Preparative HPLC Column #2, Custom Order (50 mm ID)	79644
Preparative HPLC Column #3, Custom Order (100 mm ID)	79645

From autosampler syringes to manual injection and more, Hamilton has the chromatography syringe you need. View the full portfolio of HPLC, GC and TLC syringes at [www.hamiltoncompany.com/syringes](http://www.hamiltoncompany.com/syringes).

# HPLC Guard Columns

Hamilton guard columns protect analytical, semi-prep and preparative HPLC columns. They also remove particulate contaminants and highly adsorptive compounds from samples, prolonging column life. The designs are modular and cartridge replacement is easy and tool-free.

Cartridge holders for analytical columns are available in stainless steel and PEEK materials. Semi-prep/preparative holders are available only in stainless steel. The design of the cartridge holders is optimized to reduce dead volume and prevent extra column band broadening. Cartridge holders are reusable indefinitely.

Packing Material	Analytical Starter Kits (1 holder, 2 cartridges)		Analytical Replacement Cartridges (5/pk)		Semi-prep/preparative Starter Kits (1 holder, 1 cartridge)	Semi-prep/preparative Replacement Cartridges (2/pk)
	Stainless Steel	PEEK	Stainless Steel	PEEK	Stainless Steel	Stainless Steel
PRP-1	79447	79317	79445	79318	79121	79122
PRP-3	79461	79393	79454	79395	79123	79124
PRP-C18	79685	79687	79686	79688	79689	79690
PRP-h5	79267		79268		79277	79278
PRP-X100	79448	79383	79446	79385	79125	79126
PRP-X110/ X110S	79726	79727	79728	79729		
PRP-X500		79319		79320		
PRP-X600		79361		79362		
RCX-10	79462	79378	79463	79379		
RCX-30		79371		79372		
PRP-X200	79456	79368	79449	79369	79127	79128
PRP-X400		79376		79377	79131	79132
PRP-X800	79830	79831	79832	79833		
HC-75 H <sup>+</sup>			79134		79133	79134
HC-75 CO <sup>3</sup>					79866	79865
PRP-X300	79460	79373	79453	79374	79129	79130

# Bulk Polymer Resin

Hamilton offers a range of bulk polymer resins for customers who prefer to pack their own column hardware. Chromatographic capacity and efficiency tests are conducted on both bulk and packed column resins to ensure product integrity. Samples from different lots can be purchased for method and process validation.

Please specify the quantity of resin needed in grams when ordering. Quantities from one gram to kilos are available.

Support	Support Material	Exchange Capacity	Pore Size	Particle Size						
				5 μm	7 μm	9 μm	10 μm	12-20 μm	30-50 μm	50-75 μm
PRP-C18	PS-DVB functionalized with C18	N/A	100 Å	79791						
PRP-C18	PS-DVB functionalized with C18	N/A	100 Å	79792						
PRP-C18	PS-DVB functionalized with C18	N/A	100 Å	79793						
PRP-h5	PS-DVB	N/A	300 Å	79280						
PRP-1	PS-DVB	N/A	100 Å	79578	79579		79580	79581	79582	79583
PRP-3	PS-DVB	N/A	300 Å	79701 79702						
PRP-X100	PS-DVB with Trimethylammonium Exchanger	0.19 meq/gm	100 Å	79584			79585	79586		
PRP-X500	Poly(methacrylamidopropyl Trimethylammonium chloride)	1.6 meq/gm	Superficially porous		79594			79595	79596	
PRP-X600	Poly(dimethylamidopropyl-methacrylamide)	1.6 meq/gm	Superficially porous		79597			79598	79599	
RCX-10	PS-DVB with Trimethylammonium Exchanger	0.35 meq/gm	100 Å		79703					
PRP-X200	PS-DVB Sulfonate Exchanger	35 μeq/gm	100 Å				79587	79588		
PRP-X400	PS-DVB Sulfonate Exchanger	2.5 meq/gm	N/A		79591			79592	79593	
HC-75 Ca <sup>2+</sup>	PS-DVB Sulfonate Exchanger	5.0 meq/gm	Gel-Type				79709			
HC-75 H <sup>+</sup>	PS-DVB Sulfonate Exchanger	5.0 meq/gm	Gel-Type				79711			
HC-75 Pb <sup>2+</sup>	PS-DVB Sulfonate Exchanger	5.0 meq/gm	Gel-Type				79712			



# Proper Column Care and Storage

All chromatographic columns have a finite lifetime. It is good practice to routinely monitor each column's retention characteristics and performance using appropriate analyte standards. If the column is to be stored for more than two weeks, it is necessary to use a mobile phase which will inhibit microbial growth. Solvent mobile phases containing sodium azide or high concentrations of methanol or acetonitrile are suggested.

The following precautions should be taken with Hamilton HPLC columns to achieve maximum product life:

1. Routinely monitor the column's performance.
2. Switch only between mutually miscible mobile phases.
3. Avoid the possibility of precipitation of salts in the column.
4. Use only filtered and degassed mobile phases.
5. Do not allow the column to dry out.
6. Keep the column capped with the end plugs that came with the column when not in use.
7. For prolonged storage, use a mobile phase that will inhibit bacterial and mold growth.
8. Unusually high operating pressure is an indication of a plugged inlet frit. It may be cleared by reversing flow through the column for 5–10 column volumes.
9. Using guard columns is highly recommended to remove particulate matter or impurities which may permanently bind to the polymer packing materials inside the analytical column.



### Technique Tip

When executing a wash procedure, use an appropriate flow rate based on the column's inside diameter.

Custom HPLC columns are available! From dimensions to particle size to packing materials, Hamilton can build you exactly what you need. See page 48 for more information.



### Restoring column performance

Contamination of the stationary phase from samples or eluents can cause the column performance to diminish over time. Typically, one of the following procedures will rejuvenate the performance of a column that has deteriorated.

#### Column Restoration Procedures

PRP-1, PRP-C18, PRP-3, PRP-h5	PRP-X100	PRP-X200, PRP-X300	PRP-X400	RCX-10	RCX-30
Flush with 40:40:20 (ACN:IPA:H <sub>2</sub> O)	Flush 50 mL of methanol with 1% 6 N nitric acid	Inject several times with 100 µL of 1 N nitric acid	Inject several times with 100 µL of 0.1 M potassium EDTA	Flush with 50 mL of 0.1 N sodium hydroxide	Flush with 150 mL of 0.1 N sodium hydroxide

#### Column Restoration Procedures

HC-75		HC-40	
Calcium Form	Hydrogen Form	Lead Form	Calcium Form
Flush with 1% calcium chloride at 0.1 mL/min overnight	Flush with 0.1 N sulfuric acid at 0.1 mL/min overnight	Flush with 1% lead nitrate at 0.1 mL/min overnight	Flush with 1% calcium chloride at 0.1 mL/min overnight

Reverse the column so that the flow is now entering in through what was the outlet fitting, and do not connect the outlet fitting to the detector.



### Technique Tip

Always make sure that the mobile phases are miscible and that precipitation will not occur. If necessary wash the column with a suitable intermediate solvent before changing over to the new phase. A minimum of five column volumes of the intermediate solvent should be used.



# Chromatography Syringes

Hamilton offers the most complete selection of syringes on the market for use in various applications including gas chromatography (GC) and high performance liquid chromatography (HPLC) and thin layer chromatography (TLC). Building exceptional syringes is an evolving science, which is why Hamilton is dedicated to the continuous research and development of this product line.

## HPLC and GC Autosampler Syringes

Hamilton offers a line of syringes designed to work with a wide range of the most popular autosamplers from Agilent, CTC PAL®, Spark Holland and more. Each syringe is expertly handcrafted to maximize sample integrity, process efficiency and new long-life syringe technology.


## Manual HPLC Syringes

Hamilton provides a variety of custom needle syringes that are designed for use with manual HPLC injection ports. These syringes feature an electro-polished needle that slides smoothly into the injection port, minimizing the wear on valuable seals and rotors.

## TLC Syringes

The last 19 mm (0.75 in.) toward the tip of the needle on TLC syringes is coated with PTFE. The treatment reduces the surface tension between the needle and the liquid making it ideal for reproducible sample spotting.

For more information about the Hamilton line of chromatography syringes, please visit [www.hamiltoncompany.com/syringes](http://www.hamiltoncompany.com/syringes).



# About Hamilton Company

Hamilton Company is a global enterprise with headquarters in Reno, Nevada; Franklin, Massachusetts; and Bonaduz, Switzerland and subsidiary offices throughout the world.

We are an industry leader in the design and manufacture of liquid handling, process analytics, robotics and automated storage solutions. For more than 60 years, Hamilton has been satisfying customer needs by combining quality materials with skilled workmanship to ensure the highest level of performance. Hamilton's lifelong commitment to precision and quality has earned us global ISO 9001 Certification.



Founded on the technology of analytical Microliter™ and Gastight® syringes, Hamilton has a broad offering of laboratory products including manual and semi-automated precision fluid measuring instruments, chromatography products, process sensors, laboratory electrodes, pipettes and more. Top innovations from these lines include Arc™ pH, DO and conductivity intelligent sensors, the BioLevigator™ 3D Cell Culture System, Microlab® 600 Diluters/Dispensers and the Microlab 300 Guided Pipetting System.

A pioneer in liquid handling equipment and laboratory automation technology, Hamilton Robotics is known for advancing life science and biotechnology industries through reliability, performance and flexibility. Hamilton is the industry leader in design and manufacturing with patented technologies such as Compression-induced O-Ring Expansion (CO-RE™), Total Aspiration and Dispensing Monitoring (TADM) and Anti-Droplet Control (ADC). Hamilton's platforms include Hamilton VANTAGE™, its newest vertically-integrated liquid handler, Microlab STAR, Hamilton's highest selling automated pipetting platform, and Microlab NIMBUS®, the first in its class of compact, high-speed, personalized pipetting workstations.



Hamilton Storage Technologies offers comprehensive ultra-low temperature automated sample management systems for microtube and microplate storage. Hamilton's line of biobanking and compound storage solutions, as well as consumables, are designed for a broad array of life science processes. Products include BiOS™, SAM™ and ASM™, designed for sample integrity, flexibility and reliability.

Hamilton Company is focused on blending invention and accuracy to deliver customers unparalleled products.

# HAMILTON



Your Hamilton Representative

©2013 Hamilton Company. All rights reserved.  
PAL is a registered trademark of CTC Analytics AG, Switzerland.  
All other trademarks are owned and/or registered by Hamilton Company in the U.S. and/or other countries.  
Lit. No. L80083

## HAMILTON®

Web: [www.hamiltoncompany.com](http://www.hamiltoncompany.com)

USA: 800-648-5950

Europe: +41-81-660-60-60

### Hamilton Americas & Pacific Rim

Hamilton Company Inc.  
4970 Energy Way  
Reno, Nevada 89502 USA  
Tel: +1-775-858-3000  
Fax: +1-775-856-7259  
[sales@hamiltoncompany.com](mailto:sales@hamiltoncompany.com)

### Hamilton Europe, Asia, & Africa

Hamilton Bonaduz A.G.  
Via Crusch 8  
CH-7402 Bonaduz, Switzerland  
Tel: +41-81-660-60-60  
Fax: +41-81-660-60-70  
[contact@hamilton.ch](mailto:contact@hamilton.ch)

To find a representative in your area, please visit [hamiltoncompany.com/contacts](http://hamiltoncompany.com/contacts).