

MACHEREY-NAGEL  
CHROMABOND®  
Promotion 2021



Chromatography



Modern polymeric SPE phases

- CHROMABOND® HLB  
for the extraction of medium to polar analytes  
from polar matrices
- CHROMABOND® HR-X<sub>perts</sub>  
for the extraction of hydrophobic and ionic  
analytes from polar matrices

Promotion\*  
Oct. 1<sup>st</sup> – Dec. 30<sup>th</sup>,  
2021

MACHEREY-NAGEL

[www.mn-net.com](http://www.mn-net.com)



# CHROMABOND® Promotion 2020

## Features of CHROMABOND® HLB

- Applicable for a wide range of analyte polarities
- Enhanced retention for polar compounds
- High loadability and outstanding performance
- Water wettable – even if bed runs dry, SPE can be continued
- Two particle sizes (60 µm and 30 µm) available
- The alternative to Oasis® HLB

## Features of CHROMABOND® HR-Xperts

- RP and different mixed-mode polymer phases in two particle sizes (45 µm and 85 µm)
- Broad spectrum of applications with special suitability for enrichment of pharmaceuticals from biological matrices
- High purity
- High loadability and outstanding elution properties

## Ordering information

REF	Phase	Description	Volume	Adsorbent weight	Particle size	Pack of
<b>N-vinylpyrrolidone-divinylbenzene (NVP/DVB) resin</b>						
730921	CHROMABOND® HLB	Hydrophilic-lipophilic balance	1 mL	30 mg	60 µm	30
730921P30	CHROMABOND® HLB	Hydrophilic-lipophilic balance	1 mL	30 mg	30 µm	30
730923	CHROMABOND® HLB	Hydrophilic-lipophilic balance	3 mL	60 mg	60 µm	30
730923P30	CHROMABOND® HLB	Hydrophilic-lipophilic balance	3 mL	60 mg	30 µm	30
730924	CHROMABOND® HLB	Hydrophilic-lipophilic balance	3 mL	200 mg	60 µm	30
730924P30	CHROMABOND® HLB	Hydrophilic-lipophilic balance	3 mL	200 mg	30 µm	30
730926	CHROMABOND® HLB	Hydrophilic-lipophilic balance	6 mL	200 mg	60 µm	30
<b>Polystyrene-divinylbenzene resin</b>						
730936	CHROMABOND® HR-X	Hydrophobic resin	3 mL	60 mg	85 µm	30
730936P45	CHROMABOND® HR-X	Hydrophobic resin	3 mL	60 mg	45 µm	30
730931	CHROMABOND® HR-X	Hydrophobic resin	3 mL	200 mg	85 µm	30
730931P45	CHROMABOND® HR-X	Hydrophobic resin	3 mL	200 mg	45 µm	30
730939	CHROMABOND® HR-X	Hydrophobic resin	6 mL	500 mg	85 µm	30
730956P45	CHROMABOND® HR-XC	Strong mixed-mode cation exchanger	3 mL	60 mg	45 µm	30
730952	CHROMABOND® HR-XC	Strong mixed-mode cation exchanger	3 mL	200 mg	85 µm	30
730950P45	CHROMABOND® HR-XA	Strong mixed-mode anion exchanger	3 mL	60 mg	45 µm	30
730951	CHROMABOND® HR-XA	Strong mixed-mode anion exchanger	3 mL	200 mg	85 µm	30
730735P45	CHROMABOND® HR-XCW	Weak mixed-mode cation exchanger	3 mL	60 mg	45 µm	30
730739	CHROMABOND® HR-XCW	Weak mixed-mode cation exchanger	3 mL	200 mg	85 µm	30
730747P45	CHROMABOND® HR-XAW	Weak mixed-mode anion exchanger	3 mL	60 mg	45 µm	30
730748	CHROMABOND® HR-XAW	Weak mixed-mode anion exchanger	3 mL	200 mg	85 µm	30

\* Special promotion prices only valid for the mentioned products. Offer valid from October 1<sup>st</sup> – December 30<sup>th</sup>, 2021

Registered trademarks

Oasis® Waters Corp. (USA)



ANCHEM PLUS  
Bora-Komorowskiego 56  
03-982 Warszawa  
Phone: +48 22 350 77 07  
Email: kontakt@anchemplus.pl

[www.mn-net.com](http://www.mn-net.com)

**MACHEREY-NAGEL**



MACHEREY-NAGEL GmbH & Co. KG  
Valenciennner Str. 11  
52355 Düren · Germany

DE Tel.: +49 24 21 969-0 info@mn-net.com  
CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com  
FR Tel.: +33 388 68 22 68 sales-fr@mn-net.com  
US Tel.: +1 888 321 62 24 sales-us@mn-net.com