

MACHEREY-NAGEL

Thin layer chromatography

Chromatography



Enchanting prices for small sized TLC glass plates

- Quality
- Efficiency
- Selectivity



Special offer!

Period of validity:
June 1st until
December 31st 2020

MACHEREY-NAGEL

www.mn-net.com



Since 1911

Thin layer chromatography

What is TLC?

It's a quick, easy and cost-efficient method for separating compound mixtures utilizing planar plates or sheets made of glass, aluminum or polyester coated with silica or other separating media.

What is it used for?

A common method in pharmacy, food industries, natural compound analyses or organic syntheses for control of identity, yield or composition of substances and products.

Why should customers buy MACHEREY-NAGEL TLC products?

- Comprehensive range of plate sizes, surface chemistries and backings
- Pre-coated plates ready for immediate use
- Homogeneous, smooth and well adhering layers
- Consistent high quality from batch-to-batch and from plate-to-plate

ADAMANT

Unmodified silica gel (SiOH) TLC layer on glass backing with an extra binder system for outstanding abrasion resistance. Particle size 5–17 µm.

TLC glass plates in different plate sizes for many different applications

ADAMANT TLC glass plates with UV indicator (254 nm)		
Dimensions	Glass	Price
2.5 x 7.5 cm	REF 821005	
5 x 5 cm	REF 821035	
5 x 10 cm	REF 821010	
5 x 20 cm	REF 821015	
10 x 10 cm	REF 821020	
10 x 20 cm	REF 821025	
Dimensions	Glass	Price
Suitable for TLC format up to 20 x 20 cm	REF 814019	
Suitable for TLC format up to 10 x 10 cm	REF 814018	

Good to know

The ADAMANT is the hardest TLC ready-to-use layer from MACHEREY-NAGEL.



Special offer!

Period of validity:
June 1st until
December 31st 2020



ANCHEM PLUS
Bora-Komorowskiego 56
03-982 Warszawa
Phone: +48 22 350 77 07
Email: kontakt@anchemplus.pl

www.mn-net.com

MACHEREY-NAGEL



MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6–8 · 52355 Düren · Germany

DE / International:

Tel.: +49 24 21 969-0

Fax: +49 24 21 969-199

E-mail: info@mn-net.com

CH:

Tel.: +41 62 388 55 00

Fax: +41 62 388 55 05

E-mail: sales-ch@mn-net.com

FR:

Tel.: +33 388 68 22 68

Fax: +33 388 51 76 88

E-mail: sales-fr@mn-net.com

US:

Tel.: +1 484 821 0984

Fax: +1 484 821 1272

E-mail: sales-us@mn-net.com



Since 1911